2.
- For steps 2, 3, 5, you may need some help from others. Wrap the thigh belt around the thigh. The lower limit should be above the patella of the knee.
3.
- Place the frame next to the pelvis and thigh. Locate the axis of the hinge above the greater trochanter. The top edge of the orthosis should not be above the pelvic belt. The plastic plates of the lower frame should be in the thigh area.
4.
- Pull the elastic side straps of the pelvic belt and fasten them to the pelvic belt.
5.
- Fasten the thigh straps around the thigh.

Instructions for cleaning:

Detach the padding from the stabilising splint before washing. Padding and tightening straps should be hand-washed at low temperatures (30°) with mild detergents, without softeners. Do not use bleach or strong detergents. Do not spin, do not wring, gently squeeze before drying. To dry, spread the parts of the upholstery on a towel and air-dry. Keep away from artificial sources of heat. Check the condition of the orthosis after each wash. Clean the product systematically, at least twice a month, and more often if necessary.

Composition:

PE, EVA, PA, AI, ABS, Cotton


Intended user:

The medical device is intended for both laypersons (non-professional users) and professional users.

ATTENTION!

Do not use the device if you have circulatory disorders.  
Do not use the device if it has been damaged or deformed.  
Do not attempt to repair the device yourself.  
Do not use the device of the wrong size or in a way for which it is not intended. If you experience any problems related to using the product, please get in touch with your doctor. Failure to observe proper personal hygiene and inappropriate laundry detergents can cause abrasions or other ailments associated with increased sweating and bacterial growth. Do not use the orthosis with warming gels or ointments. Any serious incident in connection with a device must be reported to the manufacturer and the competent authority of the Member State where the user and/or patient is based.

IT IS RECOMMENDED TO CONSULT A SPECIALIST BEFORE USING THE DEVICE

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HIP BRACE HIP IMMOBILISATION ORTHOSIS WITH FLEXION, EXTENSION, ABDUCTION AND ADDUCTION ADJUSTMENT RANGE

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**Instructions for use:**

- limiting adduction and abduction and controlling adduction and abduction of the leg at the hip joint;
- supporting and stabilising the joint.

**Indications:**

- hip osteoporosis;
- hip dislocation;
- after surgeries, such as arthroscopy and joint replacement.

**Contraindications:**

- hip fractures;
- conditions requiring total hip immobilisation;
- hard-to-heal wounds or skin lesions within the area of the orthosis;
- allergic reaction to the materials used.

**Device description:**

- a 3-strap fastening system,
- a flat lock with comfortable adjustment of the angular range in 15° increments, flexion from -30° to 90°, extension from -30° to 30° and immobilisation in -30°, -15°, 0°, 15°, 30° positions,
- smooth adjustment of abduction and adduction from -60° to 60°,
- telescopic adjustment of the orthosis length,
- air-permeable upholstery for user comfort,
- the orthosis is available for the left and right limb.

**Target groups of patients:**

The orthosis is designed for patients with hip dysfunction, following injury or disease, to support healing and rehabilitation by limiting adduction and abduction and controlling adduction and abduction motion of the hip;

**Disposal of the medical device:**

The device must be disposed of following the applicable local environmental regulations.

**Choosing your orthosis:**

- The first time the orthosis is fitted, it should be done upon prescription and under the supervision of qualified medical personnel. The doctor or physiotherapist should decide how much stabilisation and compression the hip needs.
- The orthosis should be selected based on the dimensions shown in the size chart.
- If the device becomes too loose or changes its position, it must be refitted starting from the first step in the section, How to put the orthosis on.

**Size chart:**

Universal size.

**How to put the orthosis on:**

- Wrap the pelvic belt around your pelvis. The top of the belt should be above the iliac crest.