3. Pull and fasten the belt around the hips on the front of the orthosis body, applying appropriate compression force. FIG. 2
4. Pull and fasten the thigh belt on the front of the orthosis body, applying appropriate compression force. FIGS. 3 and 4
5. To set the desired range of motion in the joint, release the Quick Lock (c). The lock allows immobilization of the orthosis at angles of 0°, 15°, 30°, 45°. FIG. 5
6. Press the flexion (b) and extension (a) angle lock to set the desired angle. To prevent accidental pressing of the lock, use the clamp straps. FIG. 5

Instructions for cleaning:

Before washing, detach the padding elements from the stabilizing splint. The padding elements and tension straps should be hand washed at low temperatures (30°C) with mild detergents, without using fabric softeners. Do not use bleach or strong detergents. Do not spin, squeeze gently before drying. For drying, spread the upholstery components on a towel and air dry. Keep away from artificial heat sources. Check the condition of the orthosis after each wash. Clean the product systematically, at least twice a month, or more often if necessary.

Composition:

PA, PES, PU, PE, AL




Intended user:

The medical device is intended for both laypersons (non-professional users) and professional users.

ATTENTION!

Do not use the device if you have circulatory disorders.
Do not use the device if it has been damaged or deformed.
Do not attempt to repair the device yourself.
Do not use the device of the wrong size or in a way for which it is not intended. If you experience any problems related to using the product, please get in touch with your doctor. Failure to observe proper personal hygiene and inappropriate laundry detergents can cause abrasions or other ailments associated with increased sweating and bacterial growth. Do not use the orthosis with warming gels or ointments. Any serious incident in connection with a device must be reported to the manufacturer and the competent authority of the Member State where the user and/or patient is based.

IT IS RECOMMENDED TO CONSULT A SPECIALIST BEFORE USING THE DEVICE

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EN HIP RUSH BASIC HIP BRACE WITH FLEXION ANGLE ADJUSTMENT

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**Instructions for use:**

- limitation and control of flexion and extension movements of the limb in the hip joint;
- offloading and stabilization of the joint.

Indications:

- hip osteoporosis;
- hip dislocation;
- post-operative treatment, e.g., arthroscopy, arthroplasty.

Contraindications:

- fractures within the hip joint;
- conditions requiring complete immobilization of the hip joint;
- difficult-to-heal wounds or skin lesions within the area affected by the orthosis;
- allergic reaction to the materials used.

Target groups of patients:

The orthosis is designed for patients with hip joint dysfunction, following injuries or diseases, to aid in the treatment and rehabilitation process by enabling the regulation of flexion and extension range of motion of the limb at the hip joint and compression of the soft tissues around the joint.

Disposal of the medical device:

Dispose of the medical device in accordance with local environmental protection regulations.

Choosing your orthosis:

- The orthosis should initially be fitted by and under the supervision of qualified medical personnel. The doctor or physiotherapist should determine the required level of stabilization and compression for the hip joint.
- Select the orthosis based on the dimensions presented in the size chart.
- The set includes two elastic hip belts to choose from (longer and shorter).
- If the product becomes too loose or shifts, it should be re-fitted starting from the first step of the „Fitting Instructions.”

Size chart:

Universal size for both right and left limb.

Size	Unicersal
Hip circumference (cm)	75 - 135

How to put the orthosis on:

Fitting Instructions: Initially, the orthosis is locked at a 0° position, and the length of the orthosis is 49 cm. Optimal fitting by qualified medical personnel is crucial for providing adequate stabilization of the hip joint.

- Wrap the pelvic belt around your pelvis. The top of the belt should be above Before fitting the orthosis, attach the hip and thigh belt to the rear Velcro.
- Place the orthosis on the appropriate limb, ensuring that the center of the dial is at the level of the hip joint. FIG. 1