How to put the orthosis on:

- 1. Put the orthosis on the knee so that the centre of the dials is at the level of the knee articular space.
- 2. Fasten the lower part of the orthosis and ensure that the dials and the orthosis are snug against the knee.
- 3. Fasten the upper part of the orthosis, adjusting the compression.
- 4. Fasten the additional straps: upper and lower.

Instructions for cleaning:

Detach the padding from the stabilising splint before washing. Padding and tightening straps should be hand-washed at low temperatures (30°) with mild detergents, without softeners. Do not use bleach or strong detergents. Do not spin, do not wring, gently squeeze before drying. To dry, spread the parts of the upholstery on a towel and air-dry. Keep away from artificial sources of heat. Check the condition of the orthosis after each wash. Clean the product systematically, at least twice a month, and more often if necessary.

Composition:

CR, PA, AI

Intended user:


The medical device is intended for both laypersons (non-professional users) and professional users.

ATTENTION!

It is not recommended to use the product directly on the skin. Do not use the device if you have circulatory disorders. Do not use the device if it has been damaged or deformed. Do not attempt to repair the device yourself. Do not use the device of the wrong size or in a way for which it is not intended. If you experience any problems related to using the product, please get in touch with your doctor. Failure to observe proper personal hygiene and inappropriate laundry detergents can cause abrasions or other ailments associated with increased sweating and bacterial growth. Do not use the orthosis with warming gels or ointments. Any serious incident in connection with a device must be reported to the manufacturer and the competent authority of the Member State where the user and/or patient is based.

IT IS RECOMMENDED TO CONSULT A SPECIALIST BEFORE USING THE DEVICE

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**Instructions for use:**

A universal knee orthosis with a wide range of therapeutic applications. Supports the treatment and rehabilitation of injuries and dysfunctions of the knee and surrounding areas.

Indications:

- subacute and chronic pain associated with strains and microtrauma of the ligament-muscle system,

- knee osteoarthritis and associated instability of the knee joint in the sagittal plane,
- stabilisation after surgery in the final stage of functional therapy (recovery),
- injuries and dysfunctions without immobilisation of the joint.

Contraindications:

Conditions requiring immobilisation with a plaster cast. Hard-to-heal wounds and ulcers. Significant bone deformities (inability to adapt the orthosis to the joint). Varicose veins of the lower limbs. Open wounds of the skin within the area of the orthosis.

Device description:

- an open orthosis made of neoprene,
- opening for oval popliteal fossa,
- comfortable fastening at the front of the orthosis,
- lateral stabilisation provided by polycentric stays,
- a system of transverse tightening straps.

Target groups of patients:

The orthosis is intended for patients with dysfunction of the knee, after injuries, diseases, to support treatment and rehabilitation by providing limitation of flexion and extension of the limb and lateral stabilisation of the knee.

Disposal of the medical device:

The device must be disposed of following the applicable local environmental regulations.

Choosing your orthosis:

1. The first time the orthosis is fitted, it should be done upon prescription and under the supervision of qualified medical personnel.
2. The orthosis should be selected based on the dimensions shown in the size chart.
3. If the device becomes too loose or changes its position, it must be refitted starting from the first step in the section, 'How to put the orthosis on'.

Size chart:

The orthosis is available for the left and right limb.

Size	S	M	L	XL	XXL
Thigh circumference [cm]	37–43	45–50	50–56	56–62	62–68

Thigh circumference measurement – 15 cm above the centre of the kneecap.