How to put the orthosis on:

- 1. Detach all buckles.
- 2. Place hinge axis in line with knee rotational axis. Apply the knee brace with the patient seated and knee joint flexion at 80° -90°
- 3. Press the shell against the leg and secure number 1 buckle behind the calf. Adjust the strap length to desired fit by using the hook and loop closure. Next, fasten strap number 2 around the lower calf muscles.
- 4. Fasten other straps in proper order (number 3 to 6) and in the same way as strap number 1.
- 5. Stand, and check the tightness of the straps.
- 6. Walk with the brace on to see if the strap causes any discomfort or not, or causes the pain worsening. If it does, re-adjust the strap while seated, then stand and walk again to test the pain relief.
- 7. Adjust the hinge with hex wrench for minor adjustments.

Instructions for cleaning:

Detach the padding from the stabilising splint before washing. Padding and tightening straps should be hand-washed at low temperatures (30°) with mild detergents, without softeners. Do not use bleach or strong detergents. Do not spin, do not wring, gently squeeze before drying. To dry, spread the parts of the upholstery on a towel and air-dry. Keep away from artificial sources of heat. Check the condition of the orthosis after each wash. Clean the product systematically, at least twice a month, and more often if necessary.

Composition:

Al, St, PA, PE, EVA, ABS, PU, cotton

Intended user:


The medical device is intended for both laypersons (non-professional users) and professional users.

ATTENTION!

Do not use the device if you have circulatory disorders.
Do not use the device if it has been damaged or deformed.
Do not attempt to repair the device yourself.
Do not use the device of the wrong size or in a way for which it is not intended. If you experience any problems related to using the product, please get in touch with your doctor. Failure to observe proper personal hygiene and inappropriate laundry detergents can cause abrasions or other ailments associated with increased sweating and bacterial growth. Do not use the orthosis with warming gels or ointments. Any serious incident in connection with a device must be reported to the manufacturer and the competent authority of the Member State where the user and/or patient is based.

IT IS RECOMMENDED TO CONSULT A SPECIALIST BEFORE USING THE DEVICE

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**Instructions for use:**

The orthosis is used for unicompartmental osteoarthritis of the knee, during post-operative recovery or rehabilitation periods requiring unicompartmental joint unloading, joint axis correction and stabilisation; it is also indicated for patients performing low to medium intensity activities.

Indications:

- unicompartmental osteoarthrosis of the knee/gonarthrosis,
- knee pain due to degenerative changes,
- post-traumatic reduction in the function of the musculoskeletal system stabilising the knee,
- status post ACL, PCL, MCL, LCL reconstruction,
- ACL and PCL insufficiency,
- quadriceps weakness,
- conservative treatment and stabilisation of the knee.

Contraindications:

Advanced lesions involving both compartments of the knee, conditions requiring immobilisation of the knee, wounds, ulcers, skin lesions in the area of the orthosis. Special care should be taken when using the orthosis in patients with sensory disturbances, peripheral circulation disorders, peripheral vascular conditions and diabetics.

Device description:

- an anatomical design, lightweight aluminium frame,
- smooth adjustment of valgus and varus of the orthosis frame,
- a polycentric locking system enabling dual restriction of the angular range of motion in flexion (0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°) and extension (0°, 10°, 20°, 30°, 40°) and immobilisation at 0°, 10°, 20°, 30°,
- stabilisation straps to ensure the correct positioning and effectiveness of the orthosis,
- the orthosis is available for the left and right limb.

Target groups of patients:

The orthosis is intended for patients with dysfunction of the knee, after injuries, diseases, to support treatment and rehabilitation, requiring unicompartmental unloading of the knee with the possibility of limiting flexion and extension motion of the lower limb.

Disposal of the medical device:

The device must be disposed of following the applicable local environmental regulations.

Choosing your orthosis:

- The first time the orthosis is fitted, it should be done upon prescription and under the supervision of qualified medical personnel.
- The orthosis should be selected based on the dimensions shown in the size chart.
- If the device becomes too loose or changes its position, it must be refitted starting from the first step in the section,How to put the orthosis on’.

Size chart:

The orthosis is available for the left and right limb.

Size	S	M	L	XL
Thigh circumference [cm]	39–47	47–53	53–60	60–67
Calf circumference (cm)	29–37	37–43	43–50	50–57

Thigh circumference measurement – 15 cm above the centre of the kneecap

Calf circumference measurement – 13 cm below the centre of the kneecap