Materials:

PA, PES, PU, PE, AL.

ATTENTION!

Do not use the device if you have circulatory disorders.Do not use the device if it has been damaged or deformed. Do not attempt to repair the device yourself. Do not use the device in the wrong size or in a way for which it is not intended. If you experience any problems related to using the product, please get in touch with your doctor.

ATTENTION!





Not observing the rules of personal hygiene and washing the product with inappropriate detergents may result in skin abrasion or other discomfort related to perspiration and spread of bacteria. Please do not use the product together with warming ointment or gel. After washing, rinse the product carefully - residue detergent left in the product may cause skin vexation and damage of the orthosis fabric.

ATTENTION!

Each and every serious incident related to the product should be reported to the manufacturer and competent authority of the Member State in which the user or patient is domiciled.

PLEASE SEEK MEDICAL ADVICE PRIOR TO USING THIS PRODUCT.

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Intended use:

- increasing stability of the knee joint while limiting ankle joint movement,
- limiting the range of motion of the knee joint while immobilising and stabilising the ankle joint.

Indications:

- post-knee surgery to limit movement of the ankle joint for better stability,
- knee joint injuries and neurological or orthopaedic disorders requiring stabilisation, immobilisation, or control of its range of motion combined with foot drop.

Contraindications:

- conditions requiring full limb immobilisation,
- intra-articular, tibia, fibula, femur fractures
- hard-to-heal wounds or skin lesions within the area of the orthosis;
- allergic reaction to the materials used.


Disposal of the medical device:

The device must be disposed of following the applicable local environmental regulations.

Fitting of the orthosis:

- The first time the orthosis is fitted, it should be done upon prescription and under the supervision of qualified medical personnel.
- The orthosis should be selected based on the dimensions shown in the size chart.
- If the device becomes too loose or changes its position, it must be refitted starting from the first step in the section.

Size chart:

Universal size.

Length of the orthosis from 53 to 64 cm; 49 to 66 cm (without a foot immobiliser fitted (A)).

Height of the foot immobiliser from 15 to 17.5 cm (B).

Hinge axis height from the ground from 46 to 52.5 cm (C).

Length of the orthosis with a foot immobiliser from 68 to 81.5 cm (D).

Length of the foot immobiliser 16 cm (E).

How to put the orthosis on:

- Unfold the orthosis. Place the orthosis on a flat surface.
- Then, place the orthosis under the patient's leg so that it is between the stays.
- The limb should be as shown centrally between the stays of the orthosis. Make sure the dials are correctly positioned as shown in the figure.
- Each stay has 3 grommets that can be independently moved to fit the orthosis to the limb. To move a grommet, press the button and move the grommet in the desired direction. Make sure the button is locked in the desired position, and the grommets are in the same place on both sides.
- Wrap the straps around the stays on both sides of the leg. The stays should be positioned as shown; the dial should be slightly above the knee joint. Fasten the straps in the following order: 3, 2, 4, 1. The ends of the straps can be trimmed after fitting.
- The angular range of motion (ROM) is changed by pressing the lock button and moving it to the desired range. Releasing the lock button will fix it in the correct position. Make sure the angular ranges are the same on both sides of the orthosis.
- The orthosis can be locked in 4 positions: (0, 15, 30, 45 degrees). Press and slide up the Quick Lock button to unlock and down to lock the orthosis in the chosen position. Make sure the angles on both dials are the same.
- Once the angles have been set, the dial can be secured using fastening straps. After tightening the straps, trim the end.
- Fit the foot immobiliser as follows:
 - a. tighten the two screws securing the foot immobiliser on each side (Fig. 9A) to the lower pad of the knee orthosis,
 - b. Adjust the height of the foot immobiliser. Loosen 2 screws on each side of the height adjustment (Fig. 9B), adjust to the correct height and then tighten the screws.

Instructions for cleaning:

Detach the padding from the stabilising splint before washing. Padding and tightening straps should be hand washed at low temperatures (30°) with mild detergents, without softeners. Do not use bleach or strong detergents. Do not spin, do not wring, gently squeeze before drying. To dry, spread the parts of the upholstery on a towel and air-dry. Keep away from artificial sources of heat. Check the condition of the brace orthosis each wash. Clean the product systematically, at least twice a month, and more often if necessary.

Intended user:

The medical device is intended for both laypersons (non-professional users) and professional users.