3. Put the orthosis on the patient's knee with the lower limb bent 80°–90° at the knee.
4. Fasten all buckles when the orthosis is on the limb, starting with buckle 1 around the calf. Adjust the strap length to the desired circumference using the Velcro. Fasten strap 2.
5. Fasten the thigh straps in the correct order and in the same way as strap 1.
6. The patient should stand up and check the fit of the strap circumferences.
7. The patient needs to take a few steps to check that the straps do not cause discomfort or increase pain. If so, readjust the strap in a sitting position. Repeat until the correct fit is achieved, and the patient is comfortable.
8. Make the final dial adjustment using the key provided with the orthosis.

DIAL SETTING:

1. Use a screwdriver to remove the screws and remove the dial cover.
2. Remove existing locks and insert new locks that are the desired degree of flexion and extension.
3. Screw the dial cover back on and ensure that the screws are seated and fit. Pay attention to both the inner and side hinge

Instructions for cleaning:

Detach the padding from the stabilising splint before washing. Padding and tightening straps should be hand-washed at low temperatures (30°) with mild detergents, without softeners. Do not use bleach or strong detergents. Do not spin, do not wring, gently squeeze before drying. To dry, spread the parts of the upholstery on a towel and air-dry. Keep away from artificial sources of heat. Check the condition of the orthosis after each wash. Clean the product systematically, at least twice a month, and more often if necessary.

Composition:

PP, Al, PA, PE, Cotton

Intended user:


The medical device is intended for both laypersons (non-professional users) and professional users.


ATTENTION!

Do not use the device if you have circulatory disorders.  
Do not use the device if it has been damaged or deformed.  
Do not attempt to repair the device yourself.  
Do not use the device of the wrong size or in a way for which it is not intended. If you experience any problems related to using the product, please get in touch with your doctor. Failure to observe proper personal hygiene and inappropriate laundry detergents can cause abrasions or other ailments associated with increased sweating and bacterial growth. Do not use the orthosis with warming gels or ointments. Any serious incident in connection with a device must be reported to the manufacturer and the competent authority of the Member State where the user and/or patient is based.


IT IS RECOMMENDED TO CONSULT A SPECIALIST BEFORE USING THE DEVICE

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
EN

**Instructions for use:**

The orthosis is used for unicompartmental osteoarthritis of the knee, during post-operative recovery or rehabilitation periods requiring unicompartmental joint unloading, joint axis correction and stabilisation, for use in everyday activities.

**Indications:**

- unicompartmental osteoarthritis of the knee/gonarthrosis,
- unicompartmental knee injuries requiring medial or lateral compartment unloading, i.e. articular cartilage damage, meniscus injuries, congenital dysfunction of the knee requiring medial or lateral compartment unloading and correction of the joint axis,
- knee pain due to degenerative changes,

- conservative treatment and stabilisation of the knee.

**Contraindications:**

Advanced lesions involving both compartments of the knee, conditions requiring immobilisation of the knee, wounds, ulcers, skin lesions in the area of the orthosis. Special care should be taken when using the orthosis in patients with sensory disturbances, peripheral circulation disorders, peripheral vascular conditions and diabetics.

**Device description:**

- the possibility of a perfect fit to the anatomical shapes of the limb with simultaneous application of corrective force,
- the possibility of wearing the orthosis on the skin under trousers,
- a polycentric locking system enabling dual restriction of the angular range of motion in flexion (0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°) and extension (0°, 10°, 20°, 30°, 40°) and immobilisation at 0°, 10°, 20°, 30°,
- the possibility to modify the length of the straps and the position and length of the pads, shape and adjustment of the frame,
- smooth adjustment of valgus and varus of the orthosis frame,
- stabilisation straps to ensure the correct positioning and effectiveness of the orthosis,
- the orthosis is available for the left and right limb.

**Target groups of patients:**

The orthosis is intended for patients with dysfunction of the knee, after injuries, diseases, to support treatment and rehabilitation, requiring unicompartmental unloading of the knee with the possibility of limiting flexion and extension motion of the lower limb.

**Disposal of the medical device:**

The device must be disposed of following the applicable local environmental regulations.

**Choosing your orthosis:**

The orthosis has been developed to provide maximum comfort to patients with unicompartmental gonarthrosis. The orthosis offers a perfect fit to the anatomical shapes of the limb while applying corrective force. All modifications: angular range of motion, frame shape and fit, corrective force, strap lengths, positioning and padding size should be done by qualified medical personnel.

**Size chart:**

The orthosis is available for the left and right limb.

Universal size.

**How to put the orthosis on:**

- Loosen each Velcro fastener and attach it to the strap when first fitting the orthosis.
- Unfasten all fasteners.