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Achilles heel pad

APPLICATION:

Heel pads made of EVA foam covered with nonwoven fabric, used for treatment and rehabilitation of injuries and dysfunctions of the Achilles tendon. When used with a foot and shin brace, they allow the heel to be elevated and the foot to be placed in a horse-like position, which prevents rotation. Removing the heel pads one by one allows you to gradually stretch the Achilles tendon.

CONTRAINDICATION:

Using the heel on its own, without a foot and shin brace. The heel pads cannot be used with other types of braces or with regular shoes.

PRODUCT CHARACTERISTICS:

- heel pads made of EVA foam covered with nonwoven fabric,
- five different sizes included,
- right and left versions available,
- colour: white.

APPLICATION INSTRUCTION:

1. Before attaching the heel pad, clean the inside of the walker in that area with methylated spirits and let them dry.
2. Position the heel pad in the walker.
3. Peel off the protective film of the (scratchy) dot under the heel pad. Press the heel pad firmly against the walker.
4. The heel pads can be placed on top of each other to achieve the desired heel rise. 3 heel pads equal 22° of heel raise; 2 heel pads equal 16° of heel raise; one heel pad equals 10° of heel raise.

WASHING INSTRUCTIONS:

Hands wash and wipe off. Dry in air current. Do not dry clean and bleach.

Caution: If allergic skin rashes occur, please consult your physicians.