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## INSTRUCTIONS FOR USE

### Lumbosacral brace

The Qmed Pharma lumbosacral brace effectively supports the treatment and rehabilitation process in injuries, diseases and dysfunctions of the zygapophyseal joints, intervertebral disk, muscle attachments and ligaments of the lumbar and the lumbosacral spine.

#### INDICATIONS FOR USE:

- Pain syndromes in the lumbar spine of various aetiology (discopathy, spondyloarthrosis, spondylolisthesis),
- degenerative lesions of the lumbar spine,
- hypertonia in lumbar discopathy.

#### APPLICATION:

1. Unfasten the hook-and-loop fasteners.
2. Place the belt part with the underwires against the lumbosacral spine so that the underwires expand downwards (i.e. at the top, the distances between the underwires should be smaller, while at the bottom, they should be larger).
3. Then, start from fastening the inner bottom band, lifting the underbelly (especially important in the case of large waist measurements), and then fasten the outer upper band that unburdens the spine.
4. Do not use directly on wounds or damaged skin. The skin under the belt should be clean, dry, without lotions, creams etc.

#### MAINTENANCE:

Hand-wash in water at 40°C with gentle washing agents, e.g. intended for wool. Do not let the product become very dirty. Do not chlorinate! Dry flat away from heat sources.

#### SIZE SELECTION:

In order to select a proper size, measure your waist circumference and then choose the appropriate size from the table on the packaging.

Size	S	M	L	XL	XXL
Waist (cm)	80 – 90	91 - 100	101 - 110	111 - 120	121 - 130

#### REMARKS:

Reusable product.

WHEN WEARING THE BELT, BASIC REHABILITATION EXERCISES ARE RECOMMENDED. PROLONGED BELT USE WITHOUT DOING EXERCISES MIGHT WEAKEN THE BACK MUSCLES.

THE FINAL DECISION REGARDING THE APPLICATION AND SELECTION OF EQUIPMENT IS MADE BY A SPECIALIST PHYSICIAN.