



PRODUCT
CATALOGUE
ORTHOPAEDICS

Qmed[®] **MEYRA** *medical*

We fulfill customer expectations in terms of quality, design and novelty of products and services. Cooperation with us guarantees business success to each of our partners.

Mdh Ltd offers medical solutions for orthopaedics and rehabilitation. We are inspired by our patients and their needs. Our mission is quality, design and novelty of offered products. We are proud to help people who suffer from movement disorders. Every day we improve their quality of life which gives us satisfaction and awareness that our work makes great sense. We take care about our

customers' needs – professionalism in every field is a key factor of our actions. We remain at your disposal for comprehensive consulting service. We are becoming a leader with the largest share in portfolio of goods and services purchased by medical retailers in Poland as well as an innovative supplier of non-invasive orthotic products and rehabilitation aids in Central and Eastern Europe.



Qmed is a brand of orthopaedic products owned by Mdh Ltd, a company which has existed on Polish market for more than 20 years. Mdh Ltd offers medical solution in the field of orthopaedics and rehabilitation. Company's business activities encompass manufacturing, import, export and distribution.

The company is a producer of prophylactic footwear, orthotic products and rehabilitation aids. Mdh has the largest share in portfolio of goods and services purchased by medical retailers in Poland. Owing to its consistency in strategy implementation, the company is building a position of an innovative manufacturer and a leading supplier of non-invasive orthotic and rehabilitation solutions in Europe. Mdh R&D Department is located in Lodz Special Economic Zone.

Qmed includes comprehensive orthotic solutions for people who suffer from movement disorders. The offer comprises cervical collars, trunk belts, braces for lower and upper limbs.

Qmed products effectively support treatment and rehabilitation of various motor system disorders.

All the fabrics and materials used in manufacturing process are skin-friendly and possess necessary certificates and tests. Qmed products have modern design and owing to wide range of sizes, they are very universal.

Qmed braces are useful for prophylaxis and treatment of overload and injuries of motor system dysfunctions.

Meyra is a part of Medort Group – a Polish leader of orthopaedic and rehabilitation branch. On the international market, Meyra is recognised as a legendary company in the sector of premium products. 80-years tradition, excellent products and prestige are its main asset, and now also ours.

After incorporation of Meyra to Medort Group, Mdh offer has grown to encompass all the Meyra products available on Polish market, including orthopaedic line Meyra Medical.

Apart from production of medical devices, both Meyra and Mdh companies have supported the disabled. Motto "We move people" is not only a slogan but philosophy which stands behind the business. Support and motivation of disabled people through help in gaining mobility and independence is the core of our work.

We want Meyra products to be indispensable companion of people and help them to manage little and big challenges. We consolidate the competences of different entities to find further opportunities and to implement new, innovative solutions. We are in constant contact with the disabled, doctors, clinics, therapists, caregivers and strict cooperation with qualified medical companies provide development of our products functionality.

Thanks to product trainings and workshops for retailers as well as various services, Meyra is a strong brand, which has been accompanied with experience and commitment for 80 years.



TABLE OF CONTENTS

CERVICAL COLLARS	6	●
UPPER LIMB ORTHOSES	10	●
SPINAL ORTHOSES	34	●
LOWER LIMB ORTHOSES	52	●

CERVICAL COLLARS

Cervical braces are the orthopaedic collars which support and immobilize cervical spine. Useful in treatment of muscle disorders, degenerative disease of cervical spine.

NT – BAC Collar

Schantz Type

Catalogue number: DRQA0B



PRODUCT CHARACTERISTICS:

- soft orthopaedic collar,
- light – made of polyurethane foam,
- one-piece structure,
- fastened and regulated with hook-and-loop closure,
- collar lining of anti-bacterial and anti-mycosis properties,
- anatomically profiled,
- easy to use,
- colour: white, grey.

USE:

Dysfunction of cervical area requiring relative immobilization; osteoarthritis; disc disease; other pain syndromes; follow-up treatment of bruises, twists, inflammation and strained musculoskeletal system of the cervical area.

	29,5 - 32,5	32,6 - 35,5	35,6 - 38,5	38,6 - 41,5	41,6 - 44,5	44,6 - 47,5
** 7,0	A1	A2	A3	A4	A5	A6
** 8,5	B1	B2	B3	B4	B5	B6
** 10,0	C1	C2	C3	C4	C5	C6
** 12,0	D1	D2	D3	D4	D5	D6

* neck circumference (cm) **collar height (cm)

Qmed[®] CE

NT – BAC PLUS Collar

Florida Type

Catalogue number: DRQA0C



PRODUCT CHARACTERISTICS:

- semi-rigid neck collar,
- light – made of polyurethane foam with removable, stiffening polypropylene insert,
- one-piece structure,
- fastened and regulated with hook-and-loop closure,
- collar lining of anti-bacterial properties,
- anatomically profiled,
- easy to use,
- colour: white, grey.

USE:

Dysfunction of cervical area; muscle and ligament injuries; strains and inflammation; dysfunction of vertebral joints; injury of the intervertebral disc (radical pain, referred pain); osteoarthritis; rheumatic diseases; follow-up treatment of bruises, twists, fractures.

	29,5 - 32,5	32,6 - 35,5	35,6 - 38,5	38,6 - 41,5	41,6 - 44,5	44,6 - 47,5
** 7,0	A1	A2	A3	A4	A5	A6
** 8,5	B1	B2	B3	B4	B5	B6
** 10,0	C1	C2	C3	C4	C5	C6
** 12,0	D1	D2	D3	D4	D5	D6

*neck circumference (cm) **collar height (cm)

Qmed[®] CE

Collar Stark

Rigid cervical collar

Catalogue number: MQA2D

PRODUCT DESCRIPTION:

- rigid collar of two-piece structure,
- light, made of contacted „warm“ foam,
- made of non-toxic and hypoallergenic polyurethane foam,
- fastened with hook and loop closure,
- equipped with a tracheotomy opening and additional openings for free air flow,
- X-ray, CT and MRI compatible,
- easy to use and maintain,
- colour: beige.

USE:

Immobilization of cervical area, relief for the joints, stabilization. Recommended in rehabilitation and postoperative states, for emergency use and patient's transport. Possibility of extension to the cervical-thoracic or cervical-cranial orthosis (additional construction elements).



Size	5,71 cm	8,25 cm	10,79 cm	13,33 cm
S (25-33)	S1	S2	S3	S4
M (33-40)	M1	M2	M3	M4
L (40-48)	L1	L2	L3	L4
XL (48 and more)	XL1	XL2	XL3	XL4

MEYRA
medical CE

Collar Stark Stabilizer

Catalogue number: MQA3Y

PRODUCT DESCRIPTION:

- semi-rigid construction, light and durable,
- for use with Collar Stark ,
- improves stabilisation of the cervical area in relations to the thoracic area,
- does not contain metal elements,
- X-ray, CT and MRI compatible,
- hook-and-loop closure,
- colour: beige.

USE:

Immobilization and relief of the cervical-thoracic area up to T10-T11, intended to be used with Collar Stark.

One, universal size



MEYRA
medical CE

Collar Stark Halo

Cranial ortosis Halo type

Catalogue number: MQA4L



PRODUCT DESCRIPTION:

- intended to be used with Collar Stark,
- applied on the upper part of the head it improves stabilization, restricts movement,
- does not contain metal elements,
- X ray, CT and MRI compatible,
- fastened with hook-and-loop closure,
- colour: beige.

USE:

Halo type element increases force of skull stabilisation. Used together with rigid collar makes cervical- cranial orthosis and additionally with ORRNP5 cranial-cervical-thoracic orthosis.

One, universal size.

MEYRA
medical CE

UPPER LIMB ORTHOSES

This group comprises clavicular and shoulder orthosis, elbow supports, forearm wrist orthosis, hand orthosis. The braces are great help in further treatment and convalescence in case of shoulder, elbow or wrist injury.

Acromed

Catalogue number: DRQJOB




PRODUCT CHARACTERISTICS:

- breathable fabric provides optimal comfort of use,
- system of rigid and flexible straps ensure optimal stabilization,
- two versions: left and right.

USE:

Post-operative treatment of humeral head fracture, gleoid cavity; shoulder joint dislocation; degeneration of shoulder joint and related pain syndromes, habitual dislocation of shoulder joint; follow-up treatment after reconstruction of labrumglenoidale and articular capsule; protection against shoulder joint dislocation in case of upper limb paralysis of various etiology.

Size	1	2	3	4	5
Biceps circumference at its widest point (cm)	<26	26 - 29	29 - 33	33 - 36	36 - 39



MEYRA *medical* CE

Orthopaedic vest

Catalogue number: DRQC



PRODUCT CHARACTERISTICS:

- brings arm to adducted position and into internal rotation,
- holds the elbow at 90°-120° of flexion,
- additional element to protect and relieve the hand,
- reduces the time of wearing a cast,
- fully adjustable,
- structure is similar to a classic Dessault bandage,
- made of exceptionally durable fabrics, provides maximum comfort of use,
- colour: navy blue.

USE:

Treatment; follow-up treatment of contusions and injuries of shoulder, shoulder girdle, arm and elbow joints requiring limb immobilization in adducted position and internal rotation; inflammation and degenerative changes which require immobilization and relief of upper limb; shoulder and collar-bone, or sternum and collar-bone joint dislocation.

Type	Children < 3 years	Children 3-9 years old	Children 9-14 years old	Small for adults	Medium for adults	Large for adults
Size	RO	R1	R2	R3	R4	R5
Thorax circumference (cm)	<60	<70	70 - 80	80 - 90	90 - 120	120 <

Qmed® CE

AS-BASICmedical

Arm sling

Catalogue number: MY2751

PRODUCT CHARACTERISTICS:

- padded arm sling,
- loops for thumbs and fingers,
- easy fastening by means of hook and loop fasteners with quick-lock buckles,
- usable for left and right side,
- colour: blue.

USE:

Postoperative immobilization of shoulder joint.

Size	1	2	3
Forearm length (cm)	30 – 36	36 – 42	42 – 48



MEYRA
medical CE

AS-15°medical

Arm sling with abduction cushion

Catalogue number: MY2752

PRODUCT CHARACTERISTICS:

- padded arm sling,
- anatomically formed abduction cushion,
- easy fastening by means of hook and loop fasteners with quick-lock buckles,
- loops for thumbs and fingers,
- cushion can be cleaned separately,
- usable for left and right side,
- colour: blue.

USE:

Postoperative after fractures and implantations, fracture of the humerus, luxations.

Size	1	2	3
Forearm length (cm)	30 – 36	36 – 42	42 – 48



MEYRA
medical CE



Arm abduction orthosis

Catalogue number: QSM00020 - size S
 QSM00030 - size M
 QSM00040 - size L
 QSM00050 - size XL

PRODUCT CHARACTERISTICS:

- orthosis made of polyamide knitted-fabric laminated with polyurethane foam,
- system of wedges enables shoulder immobilization at 15°, 30°, 45°, 60° of abduction,
- open construction allows to activate forearm muscles while immobilizing the shoulder,
- system of straps enables quick and easy adjustment to individual needs.

USE:

Stabilization and abduction of upper limb after injuries (sprain, dislocation) or surgeries; follow-up treatment after reconstruction of rotator cuff; immobilization and relief for upper limb in case of inflammation of shoulder joint.

Size	S	M	L	XL
Forearm length (cm)	25 - 28	28 - 30	30 - 33	33 - 36



Qmed® Silver Line CE

Arm abduction orthosis

Catalogue number: DRQJ0A

PRODUCT CHARACTERISTICS:

- light, strong frame made of steel, with a lock to stabilise the shoulder at an appropriate, therapeutic angle (15°, 30°, 45°, 60°, 75°, 90°),
- frame covered with soft lining made of breathable fabric for optimal comfort of use,
- adjustable system of stabilizing straps allows ideal fitting to any anthropometric dimensions of the patient,
- straps ensure stabilization of the abducted limb, evenly distribute its weight, prevent strain on the opposite shoulder and the cervical spine.

USE:

Treatment after shoulder dislocation; plasty or injury to the rotor cone; arthroscopy of the shoulder joint; subcapitular fracture; Bankart type injury; post-surgical stabilisation of the shoulder joint; endoprosthetics of the shoulder joint.

Size	S	M	L
Height (cm)	150 - 165	165-170	over 170



Qmed® CE

4-Arm non-woven fabric sling

Catalogue number: DRQG0A



PRODUCT CHARACTERISTICS:

- soft sling made of polypropylene non-woven fabric,
- supports elbow, forearm, hand,
- wide and long carrying strap with hook-and-loop closure ensures proper distribution of limb weight and quick sling adjustment,
- colour: black.

USE:

Dysfunctions and injuries which require support, stabilization or movement limitation of upper limb and shoulder girdle; pain syndromes, inflammation and injuries of soft tissue; support for orthopaedic cast; preventive treatment of shoulder and humeral joint dislocation in paralysis and paresis.

Size	S	M	L
Forearm length (cm)	<30	31 - 40	>41

Qmed® CE

Arm sling

Catalogue number: PHAG011



PRODUCT CHARACTERISTICS:

- universal sling,
- light and comfortable structure,
- supports the arm and enables its partial immobilization.

USE:

Support for the arm after injuries, fractures; relief after surgical intervention; helps to maintain the arm in neutral position.

Size	S	M	L	XL
Forearm length (cm)	20 - 25	25 - 30	30 - 35	35 - 40

Qmed+pharma CE

Elbow Support

Catalogue number: PHAF011

PRODUCT CHARACTERISTICS:

- multifunctional elbow support,
- helpful while practicing sport and during rehabilitation process,
- protects and stabilizes elbow joint,
- provides joint compression.

USE:

Pain relief, preventive treatment of overstains of elbow joint, help in case of joint overload.

Size	S	M	L	XL
Elbow circumference (cm)	19 - 21	22 - 24	25 - 28	29 - 32



Qmed+pharma CE

EPImedical

Elbow bandage with silicone pads

Catalogue number: MY2731


PRODUCT CHARACTERISTICS:

- knitted fabric for optimal comfort and better adjustment to body shape,
- two pads for compression distribution at the epicondyles,
- left and right version,
- colour: silver-blue.

USE:

General irritations and rheumatic diseases of the elbow, tennis and golfer elbow, distortions, arthrosis.

Size	0	1	2	3	4	5	6
Elbow circumference (cm)	17 - 19	19 - 21	21 - 23	23 - 25	25 - 27	27 - 29	29 - 31



MEYRA medical CE



Elbow Safe

Elbow support

Catalogue number: DRQF0G


PRODUCT CHARACTERISTICS:

- highly elastic, breathable fabric ensures optimal adaptation to body shape,
- additional circumferential strap with hook-and-loop closure,
- additional pads and straps relieve attachments of muscles extensor and/or flexor of carpal joint,
- ensures optimal, therapeutic compression,
- colour: silver-blue.

USE:

Follow-up treatment of elbow joint injuries; preventive treatment of overstrains of elbow joint; follow-up treatment of tennis elbow and golfer's elbow; elbow joint pain syndromes of different etiology; soft tissue disorders.

Size	1	2	3	4	5
Elbow joint circumference (measured 5 cm below the joint)	20 - 23	23 - 26	26 - 29	29 - 32	32 - 34



MEYRA *medical* CE

Elbow joint orthosis

Catalogue number: QEM061/62

PRODUCT CHARACTERISTICS:

- designed for forearm and shoulder,
- single ROM hinge to control range of motion (adjustment every 10°),
- additional carrying strap for better relief and limb support,
- length and circumference adjustment,
- two versions: left and right,
- colour: grey.

USE:

Treatment, curing and rehabilitation of injuries to joint structures which require limitation of pronation and supination of the forearm and controlled motion.

One, universal size



Qmed® Silver Line CE

Elbow Rush

Post-operative elbow brace

Catalogue number: DRQ0F



PRODUCT CHARACTERISTICS:

- single ROM hinge to control elbow joint motion-adjustment every 10°
- quick lock: 0°,15°,30°,45°
- possible elbow joint immobilization in chosen position,
- additional suspension strap to offload and support the limb,
- adjustable length and circumference,
- additional hand holder,
- two versions: left and right,
- colour: black.

USE:

Treatment, further treatment and rehabilitation of injuries which require limitation of pronation and supination of the forearm and controlled motion.

One, universal size

Brace length: 35 cm- 48 cm

Qmed® CE

Wrist support

Catalogue number: PHAF021



PRODUCT CHARACTERISTICS:

- multifunctional wrist support,
- helpful while practicing sport and during rehabilitation process,
- protects and stabilizes wrist joint,
- provides joint compression.

USE:

Pain relief, preventive treatment of overloads of wrist joint, help in case of joint overload.

Size	S	M	L	XL
Wrist circumference (cm)	15 - 16	17 - 18	19 - 20	21 - 22

Qmed+pharma CE





Prohand

Catalogue number: DRQF0F


PRODUCT CHARACTERISTICS:

- simple construction,
- made of highly elastic, breathable knitted fabric,
- hook-and-loop closure,
- universal aid for right and left wrist.

USE:

Follow-up treatment after wrist injuries; wrist inflammation; pain syndromes of different etiology; degenerative wrist disorders; wrist protection while taking up sport and professional activities after injury.

One, universal size



MEYRA
medical CE

RHIZOmedical

Thumb orthosis

Catalogue number: MY2763

PRODUCT CHARACTERISTICS:


- two hook and loop fasteners on wrist and thumb,
- shapeable metal insert at the thumb,
- fabric permeable to air,
- left and right version,
- colour: black.

USE:

Irritations on the thumb metacarpophalangeal and thumb saddle joint, skier's thumb, rhizarthritis, ligament instabilities.



Size	1	2	3
Wrist circumference (cm)	12,5 – 16	16 – 20	20 – 23



MEYRA
medical CE



MANUmeyra Comfort

Soft wrist bandage

Catalogue number: MY2765

PRODUCT CHARACTERISTICS:


- adjustable to body shape,
- soft, skin-friendly,
- bi-elastic knitted fabric without constrictions at the edges,
- hook and loop fastener,
- anatomically shaped metal insert for optimal stabilization,
- elastic strap for optimal stabilization,
- two versions available: left and right.

USE:

Soft-tissue irritations of the wrist, immobilisation after distortions and in case of tendovaginitis, carpal tunnel syndrome and rheumatoid arthritis.



Size	1	2	3	4	5	6
Wrist circumference (cm)	14 – 15	15 – 16	16 – 17	17 – 18	18 – 19	19 – 20



MEYRA
medical CE

MANUmedical

Wrist orthosis

Catalogue number: MY2762

PRODUCT CHARACTERISTICS:


- stable orthosis for wrist positioning,
- adjustable by means of three hook and loop straps,
- stable, shapeable aluminium insert,
- fabric permeable to air,
- left and right version,
- colour: black.

USE:

Injuries of the thumb saddle joint, irritations in the wrist area, arthrosis and distortions, rheumatoid arthritis, instabilities.



Size	1	2	3	4
Wrist circumference (cm)	13 – 16	16 – 19	19 – 22	22 – 25



MEYRA
medical CE





VEOLIA

EVENTUS
zdrowie i aktywność



RUGBY
BUDOWLANSKI



Manu Universal

Wrist orthosis with a thumb hold

Catalogue number: DRQF3C



PRODUCT CHARACTERISTICS:

- light, open structure made of foam, nylon and polyester,
- anatomically profiled aluminium stay,
- additional stay with thumb spica to stabilize the thumb,
- anatomical shape allows to immobilize the thumb in functional position,
- three tightening straps enable proper adjustment,
- soft fabrics and trim guarantees optimal comfort of use,
- universal for left and right hand,
- colour: black.

USE:

Treatment of wrist, thumb and distal forearm, chronic instability of wrist and thumb with pain syndrome, follow-up treatment after surgery or conservative treatment, wrist and thumb protection while practicing sport or performing jobs which may lead to injury or overload.

One, universal size.

Qmed® CE

Hand and forearm ortosis

Catalogue number: Q41301



PRODUCT CHARACTERISTICS:

- light, made of foam, nylon, polyester and lycra,
- two splints for additional reinforcement,
- enables unhampered thumb movement,
- ensures wrist joint immobilization in functional position,
- removable stay of anatomical shape,
- universal model for left and right hand,
- colour: black.

USE:

Degenerative changes and overload of wrist joint, pain and inflammation, wrist protection while practicing sport or performing jobs which may lead to injury.

Size	S	M	L	XL
Wrist circumference (cm)	13-15	15-18	18-20	20-22

Qmed® Silver Line CE

Versa Fit

Hand and forearm orthosis

Catalogue number: MQC2E

28



PRODUCT CHARACTERISTICS:

- light, open construction made of foam, nylon and polyester fabric,
- anatomically moulded, removable, aluminium stay reinforced with two splints,
- enables wrist joint immobilization in functional position,
- three tightening straps ensure proper adjustment,
- soft fabrics provide optimal comfort of use,
- two versions available: left and right,
- colour: black.

USE:

Treatment of wrist, distal part of forearm, chronic wrist instability with pain syndrome, follow-up treatment, wrist protection while practicing sport or performing work which may lead to injuries or/in overload.

Size	XS	STANDARD	XL
Wrist circumference (cm)	<9	9-25	>22

MEYRA
medical CE

Thumb Versa Fit

Hand and forearm orthosis
with a thumb hold

Catalogue number: MQC2C



PRODUCT CHARACTERISTICS:

- light, open construction made of foam, nylon and polyester fabric,
- anatomically moulded, removable, aluminium stay reinforced with two splints,
- specially moulded aluminium stay enables thumb stabilization,
- ensures wrist joint immobilization in functional position,
- three tightening straps ensure proper adjustment,
- soft fabrics provide optimal comfort of use,
- two versions available: left and right,
- colour: black.

USE:

Treatment of wrist, thumb and distal part of forearm, chronic wrist and thumb instability with pain syndrome, follow-up treatment, wrist and thumb protection while practicing sport or performing work which may lead to injuries or/in overload.

Size	XS	STANDARD	XL
Wrist circumference (cm)	<9	9-25	>22

MEYRA
medical CE

MANUmedical Plus

Wrist orthosis with a thumb hold

Catalogue number: MY2761




PRODUCT CHARACTERISTICS:

- stable orthosis for wrist positioning,
- adjustable by means of three hook and loop straps,
- with stable, shapeable aluminium insert,
- fabric permeable to air,
- left and right version,
- colour: black.

USE:

Injuries of the thumb saddle joint, irritations in the wrist area, arthrosis and distortions, rheumatoid arthritis, instabilities.

Size	1	2	3	4
Wrist circumference (cm)	13 – 16	16 – 19	19 – 22	22 – 25



MEYRA
medical CE

Uni Hand

Splint for arm and forearm

Catalogue number: DRQF3D



PRODUCT CHARACTERISTICS:

- light, open structure padded with soft, hypoallergenic knitwear which guarantees comfort of using,
- enables wrist immobilization in functional position,
- removable, anatomically profiled stay,
- fastened with hook and loop fastener,
- universal version for right and left hand,
- colour: black.

USE:

Treatment of wrist, metacarpus, fingers which require immobilization, stabilization and offload, surgery, tendonitis, carpal tunnel syndrome, peripheral nerves injuries and central nervous system (with limb paralysis), chronic polyarthritis, prevention of contractures.

Size	S	M	L
Wrist circumference (cm)	12,5 - 15,5	15,5 - 17,5	17,5 - 20,5
Hand width (cm)	8	9	10

Qmed[®] CE

Splint for palm and forearm with thumb spica

30

Catalogue number: Q42006



PRODUCT CHARACTERISTICS:

- wrist and forearm splint with thumb spica,
- fixed with wide straps with hook-and-loop closure,
- left and right version.

USE:

Diseases of wrist, hand and fingers which require support and stabilization of dressed area in particular corrective positions, trigger fingers, carpal tunnel syndrome, rheumatoid arthritis, prevention of contractures, peripheral nerve injuries, surgeries.

Size based on hand measurement at the height of metacarpophalangeal joints.

Size	XS	S	M	L
Hand width (cm)	6 - 7	7 - 8	8 - 9	9 - 10

Qmed® Silver Line CE

WHOSP - F

Wrist and forearm splint

Catalogue number: DRQF0C



PRODUCT CHARACTERISTICS:

- wrist and forearm splint,
- adjustable to patient's needs and progress, light structure made of steel and aluminium covered with PVE foam,
- removable, soft padding with fingers separator,
- hook-and-loop closure,
- universal version.

USE:

Wrist, hand and finger affections requiring support and stabilization of the affected area in specified corrective and/or pain reducing settings, sinew inflammation, carpal tunnel syndrome, rheumatic joint inflammation, contractures, preventive treatment, circumferential nerves injuries, surgery.

Size based on hand measurement at the height of metacarpophalangeal joints.

Size	S	M	L
Hand width (cm)	6,5 - 8,0	8,0 - 9,5	>9,5

Qmed® CE

WHOSP - FT

Splint for palm and forearm with thumb spica

Catalogue number: DRQF0D



PRODUCT CHARACTERISTICS:

- wrist and forearm splint with thumb spica,
- adjustable to patient's needs and progress, light structure made of steel and aluminium covered with PVE foam,
- removable, soft padding with fingers separator,
- left and right version.

USE:

Wrist, hand and fingers affections requiring support and stabilization of the supplied area in specified corrective and/or pain reducing settings, sinew inflammation, carpal tunnel syndrome, rheumatic joint inflammation, contractures, preventive treatment, circumferential nerves injuries, surgery.

Size based on hand measurement at the height of metacarpophalangeal joints.

Size	S	M	L
Hand width (cm)	6,5 - 8,0	8,0 - 9,5	>9,5

Qmed® CE

Stack apparatus

Catalogue number: DRQF3A



PRODUCT CHARACTERISTICS:

- light, hygienic finger splint,
- highly effective,
- colour: beige.

USE:

Immobilization of distal interphalangeal joints, finger extensor injury, avulsion fracture, fracture of distal phalanx.

Size	2	3	4	5	5,5	6	7
Circumference of distal interphalangeal joint (cm)	5,5-6,0	6,0-6,5	6,5-7,0	7,5-8,0	8,5-9,0	9,5-10,0	10,5-11,0

Qmed® CE

SPINAL ORTHOSES

Braces which support cervical and lumbosacral spine. The group comprises rigid or semi-rigid corsets useful in treatment of acute or chronic pain.

Abdominal binder

Catalogue number: DRQH6



PRODUCT CHARACTERISTICS:

- made of flexible, breathable fabric,
- wide hook and loop fastener,
- perfectly adjusts to body shape,
- provides optimal tissue compression,
- colour: white.

USE:

Protection and stabilization after surgery, in postpartum period.

Size	S	M	L	XL	XXL
Waist circumference (cm)	56 - 71	69 - 84	81 - 96	94 - 109	107 - 122

Height: 26 cm

Qmed® life CE

Support stoma waistband

Catalogue number: DRQH0T



PRODUCT CHARACTERISTICS:

- made of elastic fabrics,
- equipped with ostomy opening,
- fastened with hook and loop fastener,
- wide range of sizes,
- colour: white.

USE:

Designed for patients with hernia, reduces risk of hernia in case of intensive fistula (especially in post-operative period). Strengthens stomach muscles during sport activities.

Size	XS	S	M	L	XL	XXL	XXXL
Abdominal circumference (cm)	<80	81-90	91-100	101-110	111-120	121-140	>140

Height: 24 cm

Qmed® life CE





Lumbosacral corset

Catalogue number: PHAH011

PRODUCT CHARACTERISTICS:

- effectively supports rehabilitation process and treatment,
- three baleens in the back ensure external stabilization,
- two tightening straps to adjust compression,
- colour: beige.

USE:

Injuries and dysfunctions of following: intervertebral joints, intervertebral disc, muscle attachments and ligaments of lumbar spine; pain syndromes of various etiology (discopathy, spondylosis, spondylolisthesis); degenerative changes of lumbar spine, increased muscular tension in case of discopathy.

Size	S	M	L	XL	XXL
Waist circumference (cm)	80 - 90	91 - 100	101 - 110	111 - 120	121 - 130

Height: 24 cm



Qmed+pharma CE

LS Support

Lumbosacral orthosis

Catalogue number: QBE010

PRODUCT CHARACTERISTICS:

- made of flexible strap,
- two tightening straps to regulate dorsal compression,
- five metal stays enable passive stabilization.

USE:

Nonoperativetreatment of lumbosacral area dysfunctions of different etiology; stabilization of lumbar area in case of segmental instability after surgeries; follow-up treatment after neurosurgical/orthopaedic interventions (functional treatment); prophylaxis against injuries, overload; congenital defects of lumbar area (with or without pain syndromes); degenerative spinal disorders with pain syndrome.

Size	S	M	L	XL	XXL
Waist circumference (cm)	56 - 71	69 - 84	81 - 96	94 - 109	107 - 122

Height: 24 cm



Qmed® Silver Line CE

Spincare Flex III

Lumbosacral orthosis

Catalogue number: DRQH0B1FG



PRODUCT CHARACTERISTICS:

- made of flexible, tensile straps,
- four semi-elastic, removable stays provide passive stabilization,
- breathable fabric enables warmth and moisture to be carried away,
- colour: beige.

USE:

Non-invasive treatment in dysfunctions of the lumbosacral segments of the various etiology and the degree of the intensity, stabilization of the lumbar region in case of the segment instability after surgeries, curing after neurological and orthopaedic surgeries, functional treatment, prevention of another injuries, overstrains, inborn defects of the lumbar segment of the lumbar (with or without pain syndromes), chronic pain ailments in the course of the degenerative disease.

Size	S	M	L	XL	XXL
Waist circumference (cm)	56 - 71	69 - 84	81 - 96	94 - 109	107 - 122

Height: 27 cm

Qmed® CE

Spinactive

Lumbosacral orthosis

Catalogue number: DRQH0A



PRODUCT CHARACTERISTICS:

- made of flexible, tensile straps,
- four semi-elastic stays provide passive stabilization,
- breathable material enables warmth and moisture to be carried away,
- additional removable straps allow to adjust dorsal compression and the direction of dynamic posture correction,
- dynamic correction of posture,
- colour: beige.

USE:

Intense lumbar pain, disc disease, sciatica, painful degeneration of spine, straining syndromes, stabilization of lumbar area after surgery during functional treatment, disorders of abdominal muscle tension (surgical intervention, risk of ruptures).

Size	1	2	3	4	5	6
Waist circumference (cm)	52 - 62	63 - 74	75 - 89	90 - 105	106 - 120	121 - 140

Height: 25,5 cm

Qmed® CE

Spinactive - LS II

Lumbosacral orthosis

Catalogue number: DRQHOK

PRODUCT CHARACTERISTICS:

- made of flexible, tensile straps,
- four semi-elastic stays provide passive stabilization,
- breathable fabric enables warmth and moisture to be carried away,
- removable straps allow to adjust dorsal compression,
- colour: grey.

USE:

Intense lumbar pain, disc disease, sciatica, painful degeneration of spine, straining syndromes, stabilisation of lumbar area after surgery during functional treatment, disorders of abdominal muscle tension (surgical intervention, risk of ruptures).

Size	1	2	3	4	5	6
Waist circumference (cm)	52 - 62	63 - 74	75 - 89	90 - 105	106 - 120	121 - 140

Height: 25,5 cm



Qmed® CE

Spincare Frame

Lumbosacral orthosis

Catalogue number: DRQH0I

PRODUCT CHARACTERISTICS:

- made of elastic fabric,
- enables warmth and moisture to be carried away,
- four elastic stays provide passive stabilization,
- special belt loops facilitate fastening,
- colour: grey.

USE:

Intense pain of lumbosacral spine, disc disease, stabilization after surgery and during functional treatment, instability of lumbosacral area.

Size	1	2	3	4	5
Waist circumference (cm)	52 - 62	63 - 74	75 - 89	90 - 105	106 - 120

Height: 26 cm



Qmed® CE

LUMBOmedical

Lumbar bandage with sacrum pad

Catalogue number: MY2741




PRODUCT CHARACTERISTICS:

- knitted fabric for optimal comfort and better adjustment to body shape,
- sacrum pad for the relief and massage of the surrounding musculature,
- steady and convenient hook and loop fastener with loops,
- colour: silver-blue.

USE:

Lumbagos, irritations of the iliosacral joints, muscular debility, attrition, backache, ligamentosis.

Size	1	2	3	4	5	6
Waist circumference (cm)	70 – 80	80 – 90	90 – 100	100 – 110	110 – 120	120 – 130



MEYRA
medical CE

LUMBOmedical Comfort

Lumbar bandage with sacrum pad

Catalogue number: MY2742




PRODUCT CHARACTERISTICS:

- knitted fabric for optimal comfort and better adjustment to body shape,
- sacrum pad for relief and massage of the surrounding musculature,
- steady and convenient hook and loop fastener with loops,
- colour: silver-blue.

USE:

Lumbagos, irritations of the iliosacral joints, muscular debility, attrition, backache, ligamentosis.

Size	1	2	3	4	5	6
Waist circumference (cm)	70 – 80	80 – 90	90 – 100	100 – 110	110 – 120	120 – 130



MEYRA
medical CE

LUMBOmedical PRO

Lumbar bandage with 4 stabilization stays

Catalogue number: MY2743


PRODUCT CHARACTERISTICS:

- knitted to fit the anatomical forms and provide high wearability,
- 4 stays provide more stability in the lumbar area,
- steady and convenient hook and loop fastener with loops,
- colour: silver-blue.

USE:

Ligamentosis, symphysiolysis, lumbagos, irritations of the intervertebral joints, backache, lordosis.

Size	1	2	3	4	5	6
Waist circumference (cm)	70 – 80	80 – 90	90 – 100	100 – 110	110 – 120	120 – 130



MEYRA
medical CE

LUMBOmedical Light

Lumbar support bandage
with 4 stabilisation stays

Catalogue number: MY2744


PRODUCT CHARACTERISTICS:

- fabric elastic and permeable to air,
- fixation by means of hook and loop fasteners and elastic straps,
- four reinforcing stays in the lumbar area,
- adjustable elastic stays,
- hook and loop fastener with loops.

USE:

Lumbagos, lordosis, symphysiolysis, attrition, general backache in the lumbar area.

Size	1	2	3	4	5	6
Waist circumference (cm)	64 – 72	72 – 80	80 – 90	90 – 100	100 – 110	110 – 120



MEYRA
medical CE

Premium S

Lumbosacral orthosis
Catalogue number: DRQH0R



PRODUCT CHARACTERISTICS:

- made of soft, flexible, air permeable and hypoallergenic fabric,
- two semi-rigid stays, two reinforcing elements in the dorsal part and two in the abdominal part ensure external stabilization of the sacro-lumbar section,
- additional tightening straps allow adjustment of the dorsal compression according to the level of injury and control of the direction of the dynamic correction of posture.

USE:

Non-invasive treatment of the movement segment of different etiology and different degree of intensity, stabilisation of the lumbar segment in the cases of segmentary instability following surgeries, follow-up treatment after neurological and orthopedic surgeries – functional treatment, prophylaxis against recurrent injuries, occupational and sporting excessive loads, inborn faults of the lumbar segment of the spine (with and without pain syndromes), chronic pain in degenerative diseases.

Size	S	M	L	XL	XXL
Waist circumference (cm)	85 - 95	95 - 105	105 - 115	115 - 125	125 - 135

Height: 29 cm

Qmed® CE

Spincare PIN

Lumbosacral orthosis
Catalogue number: QBE230



PRODUCT CHARACTERISTICS:

- main body made of a soft and breathable double knitted fabric,
- intuitive lacing system to adjust compression to individual needs,
- polyamide lace to manage the tension up to 16 kg,
- wide range of sizes,
- colour: black.

USE:

Nonoperative treatment of spinal segment dysfunctions of different etiology; stabilization of lumbar area in case of segmental instability after surgeries; follow-up treatment after neurosurgical/orthopedic interventions; prophylaxis against injuries, overload; congenital defects of lumbar area (with or without pain syndromes); degenerative spinal disorders with pain syndrome.

Size	S	M	L	XL	XXL	XXXL
Waist circumference (cm)	56-71	69-84	81-96	94-109	107-122	120-135

Height: 24 cm

Qmed® Silver Line CE

Spincare Flex

Thoracic – Lumbosacral orthosis

Catalogue number: DRQH0B1AA

PRODUCT CHARACTERISTICS:

- dynamic control of movement range,
- high stabilization,
- five semi-elastic asymmetric stays provide passive stabilization,
- optimal anatomical adjustment,
- perfect air and moisture permeability,
- adjustable height and stabilization strength of TLS section,
- highly comfortable,
- colour: beige.

USE:

Intense spinal pain, disc disease, painful degeneration of spine, stabilization after surgery during functional treatment, spinal stenosis.

Size	1	2	3	4	5	6
Waist circumference (cm)	52 - 62	63 - 74	75 - 89	90 - 105	106 - 120	121 - 140

Height: 31 cm



Qmed® CE

Spincare Flex II

Thoracic – Lumbosacral orthosis

Catalogue number: DRQH0B1LA

PRODUCT CHARACTERISTICS:

- made of flexible straps,
- knitwear in the dorsal part and two oblique straps to support dorsal muscles,
- additional straps to adjust dorsal compression,
- four semi-elastic, asymmetric stays provide passive stabilization,
- colour: black.

USE:

Nonoperative treatment of spinal segment dysfunctions of different etiology, stabilization of lumbar spine in case of segmental instability after surgeries, follow-up treatment after neurosurgical/orthopedic interventions, prophylaxis against injuries, overload, congenital defects of lumbar spine (with or without pain syndromes), degenerative spinal disorders with pain syndrome.

Size	1	2	3	4	5	6
Waist circumference (cm)	52 - 62	63 - 74	75 - 89	90 - 105	106 - 120	121 - 140

Height: 31 cm



Qmed® CE

Corset Premium

Thoracic-Lumbosacral orthosis

Catalogue number: DRQHOU



PRODUCT CHARACTERISTICS:

- made of soft, flexible, air permeable and hypoallergenic fabric,
- two semi-rigid plastic stays and two baleens in the back part and two baleens in the abdominal part (placed in the tightening belts) ensure external stabilization of the sacro-lumbar section,
- additional tightening straps to adjust dorsal compression,
- colour: grey.

USE:

Non-invasive treatment of the movement segment of different etiology and different degree of intensity; stabilisation of the thoracic-lumbar-sacral segment in the cases of segmentary instability following surgeries; follow-up treatment after neurological and orthopedic surgeries – functional treatment; prophylaxis against recurrent injuries, occupational and sporting excessive loads; inborn faults of the thoracic-lumbar-sacral segment of the spine (with and without pain syndromes); chronic pain in degenerative diseases.

Size	S	M	L	XL	XXL
Waist circumference (cm)	85 – 95	95 – 105	105 – 115	115 – 125	125 – 135
Height (cm)	33	33	35	35	35

Qmed[®] CE

Spincare Rigid II

Thoracic – Lumbosacral orthosis

Catalogue number: DRQH0J



PRODUCT CHARACTERISTICS:

- made of flexible, tensile knitted fabric,
- the material carries moisture away from the skin surface and in the place of orthosis adherence, allows to dissipate the excessive heat,
- four metal, profiled stays provide passive stabilization,
- additional elements for lateral stabilization,
- two additional flexible bands allow to adjust direction and compression strength.

USE:

Prevention and treatment of pain of thoracic-lumbosacral area (distribution of compression zones on the dorsal spine), degenerative spinal disorders with pain syndrome, stabilization after surgery during functional treatment, support for the abdominal wall while lowering the bowels, the structure of the orthosis allows for adjustment of the strength and height of stabilisation to the individual needs, control of TLS section movement.

Size	1	2	3	4	5
Waist circumference (cm)	52 – 62	63 – 74	75 – 89	90 – 105	106 – 120

Height: 33 cm (size 1-3)
35 cm (size 4-5)

Qmed[®] CE

Spincare High

Thoracic – Lumbosacral orthosis

Catalogue number: DRQHON

44



PRODUCT CHARACTERISTICS:

- made of elastic, flexible straps,
- high metal stays and reinforcing part,
- optimal anatomical adjustment,
- adjustable stabilization strength,
- additional straps for stabilization and posture correction,
- very good air and moisture permeability,
- highly comfortable.

USE:

Lower back pain, disc diseases – degenerative spinal disorders with pain syndromes, stenosis of vertebral canal, post-surgery stabilization during functional treatment, thoracic kyphosis, Scheuermann's disease, osteoporosis, posture defects, roundback.

Size	S	M	L	XL	XXL	XXXL
Waist circumference (cm)	80 - 90	90 - 100	100 - 110	110 - 120	120 - 130	130 - 140

Height of lumbosacral brace: 30 cm

Height of lumbosacral brace without upright keeper: 54 cm

Qmed[®] CE

Spincare High Plus

Thoracic-Lumbosacral orthosis

Catalogue number: DRQHOP



PRODUCT CHARACTERISTICS:

- made of elastic, tensile straps,
- consists of two parts: lumbosacral brace and upright keeper (hook-and-loop closure),
- two metal, profiled stays, two semi-elastic stays and two reinforcing parts placed in the lumbosacral part,
- two easy to profile stays placed in the upright keeper,
- soft armpit lining prevents abrasions,
- two bands allow to adjust compression,
- two ways of using: only the lumbosacral brace or with the upright keeper.

USE:

Intense and chronic lumbar pains, slipped disc, disc diseases, spinal stenosis, stabilization after surgery and during functional treatment, thoracic kyphosis, Scheuermann's disease, osteoporosis, weakness of abdomen muscles, faulty posture, ankylosing spondylitis.

Size	S	M	L	XL	XXL
Waist circumference (cm)	80 - 90	90 - 100	100 - 110	110 - 120	120 - 130

Height of lumbosacral brace: 31 cm

Height of lumbosacral brace without upright keeper: 54 cm

Qmed[®] CE

Premium TLS

Thoracic-Lumbosacral orthosis

Catalogue number: DRQH05



PRODUCT CHARACTERISTICS:

- made of flexible, laminated foam covered with knit-mesh out of soft polyamide Meryl-type fibre on one side and cotton knit on the other,
- four elastic dorsal stays provide best fit, motion limitation and relief,
- two additional stays ensure external stabilization of the lumbosacral part,
- four reinforcing elements in the abdominal part provide optimal stabilization,
- additional straps to adjust compression according to patient's needs,
- shoulder straps help to correct and keep proper posture.

USE:

Multilevel disc disease; spine instability; chronic pain; degenerative diseases; TLS section pain syndromes of different etiology; stabilization after surgeries; follow-up treatment.

Size	S	M	L	XL	XXL
Waist circumference (cm)	85 - 95	95 - 105	105 - 115	115 - 125	125 - 135

Height: 51 cm





Hyperextension brace Jewett type

Catalogue number: DRQB0D



PRODUCT CHARACTERISTICS:

- light, aluminium construction which enables optimal adjustment to individual needs,
- soft lining provides optimal comfort of use,
- width/height adjustable upper part of the corset,
- easy and comfortable fastening which adjusts dorsal compression,
- available also with V stabilizer.

USE:

Stabilization of thoracic-lumbosacral area (T6 - L - S) in violent, chronic and inoperable conditions, reclinacion of vertebral trunks, stabilization of T6 - L - S area after neurosurgery, oncologic surgery, and others.

Size	S	M	L
Distance between sternum manubrium to pubic symphysis (cm)	37 - 44	41 - 48	45 - 52

Qmed[®] CE

Hyperextension brace Jewett type

Movable pelvic band
Catalogue number: DRQB0D



PRODUCT CHARACTERISTICS:

- light, aluminium construction which enables optimal adjustment to individual needs,
- soft lining provides optimal comfort of use,
- width/height adjustable upper part of the corset,
- easy and comfortable fastening which adjusts dorsal compression,
- fitted for sitting position.

USE:

Stabilization of thoracic-lumbosacral area (T6 - L - S) in violent, chronic and inoperable conditions, reclinacion of vertebral trunks, stabilization of T6 - L - S area after neurosurgery, oncologic surgery, and others.

Size	S	M	L
Distance between sternum manubrium to pubic symphysis (cm)	37 - 44	41 - 48	45 - 52

Qmed[®] CE

Clavicle support

Catalogue number: QSM00091

PRODUCT CHARACTERISTICS:

- made of knitted fabric laminated with foam,
- wide range of adjustment possibilities allow to correct the body posture in patients at different ages.

USE:

Collar bone injuries which do not require surgical and orthopedic intervention, postural defects which do not require surgical intervention, pain syndromes of thoracic spine, thoracic hyperkyphosis.

One, universal size



Qmed® Silver Line CE

Spincare Extensor

Upright keeper

Catalogue number: DRQH0L

PRODUCT CHARACTERISTICS:

- made of elastic straps,
- four semi-flexible stays for stabilization,
- breathable fabric of very good air and moisture permeability,
- equipped with inserts for armpit protection.

USE:

Collar-bone injury which does not require surgical or orthopaedic intervention, posture defects which do not require surgical intervention, pain syndromes of thoracic section of the spine, kyphosis.

Size	1	2	3	4	5	6
Chest circumference (cm)	40 - 50	50 - 60	60 - 70	70 - 80	80 - 90	90 - 100



Qmed® CE



LOWER LIMB ORTHOSES

Braces which support, immobilize and protect knee and ankle joint. The group comprises simple compression braces and post-operative braces which control flexion/extension movement.

Patella tendon strap

Catalogue number: QKN00091



PRODUCT CHARACTERISTICS:

- soft knitted fabric and silicon pad to compress the patella tendon,
- regulated band allows to choose optimal circumference of the strap and proper compression.

USE:

Patellofemoral syndrome, inflammation of patella ligament, Osgood-Schlatter disease, jumper's knee.

One, universal size

51

Qmed[®] Silver Line CE

Knee support

Catalogue number: PHA1021



PRODUCT CHARACTERISTICS:

- multifunctional knee support,
- protects and supports knee joint,
- ideal for practicing sport,
- helpful in rehabilitation process.

USE:

Strain injury, reinforcement of knee joint, protection against injuries.

Size	S	M	L	XL
Knee joint circumference (cm)	27 - 33	30 - 38	39 - 45	46 - 49

Qmed+pharma CE





Knee orthosis

Catalogue number: DRQI5B

54



PRODUCT CHARACTERISTICS:

- consists of tutor made of multidirectionally flexible, breathable knitted fabrics,
- silicone ring around the patella reduces its lateral mobility but enables full range of knee movement,
- one, universal model for left and right limb.

USE

Prevention of strain injuries, pain syndromes of various etiology, inflammation (also with swelling), stabilization and relief of muscle and ligament structures in case of degenerative changes, knee protection during sport and professional activities, relief and stabilization of patella.

Size	S	M	L	XL
Knee circumference measured at the patella midpoint level (cm)	33 – 36	36 – 38	38 – 41	41 – 43

Qmed® CE

GENUmeyra

Patella bandage with silicone ring

Catalogue number: MY2711




PRODUCT CHARACTERISTICS:

- knitted fabric for optimal comfort and better adjustment to body shape,
- patella silicone ring for decompression,
- lateral coil springs for more safety,
- elastically knitted edges without constrictions,
- left and right version,
- colour: silver-blue.

USE:

Ligament injuries and insufficiencies, postoperative and traumatic irritations, joint effusions and swellings, slight instabilities, gonarthrosis.



MEYRA medical CE

Size	0	1	2	3	4	5	6	7	3 XL	4 XL	5 XL	6 XL	7 XL
Shin circumference (cm)	25 – 28	28 – 31	31 – 34	34 – 37	37 – 40	40 – 43	43 – 46	46 – 49	34 – 37	37 – 40	40 – 43	43 – 46	46 – 49
Thigh circumference (cm)	35 – 38	38 – 41	41 – 44	44 – 47	47 – 50	50 – 53	53 – 56	56 – 59	47 – 50	50 – 53	53 – 56	56 – 59	59 – 62

Motive

Open knee brace

Catalogue number: DRQ10C



PRODUCT CHARACTERISTICS:

- open knee brace made of elastic, breathable spacer fabric AirX®,
- equipped with patellar opening,
- covered popliteal fossa,
- stays with polycentric hinge provide lateral stabilization,
- transverse straps provide additional stabilization,
- colour: black.

USE:

Repetitive, subacute pain related to strain injury and micro-injury of muscles and ligaments, gonarthrosis and knee joint instability, stabilization after injury during the final phase of functional treatment, injury and dysfunction which do not require knee immobilization (lateral ligament injuries).

Size	1	2	3	4	5
Thigh circumference* (cm)	36 - 39	40 - 43	44 - 48	49 - 54	55 - 60

* In order to choose a proper size, measure the thigh circumference 15 cm above the center of patella

Length: 32,5 cm

Qmed® CE

Reflex

Knee ROM brace

Catalogue number: Wrap-around version: MQB3K

Closed version: MQB3M



PRODUCT CHARACTERISTICS:

- made of breathable spacer fabric which carries humidity away and provides optimal comfort of use,
- lateral splints may be adjusted to anatomical shape of the limb which allows to control range of motion and improves comfort of use,
- removable straps may be applied on right and left side and their marking facilitates proper fastening,
- delicate, crease-resistant fabric of specially chosen tension protects against irritation and compression of popliteal fossa,
- closed version equipped with patella stabilizer ring which location and length can be adjusted according to patient's need,
- colour: black.

USE:

Mild knee instability of various etiology.

Size	XS	S	M	L	XL	XXL	XXXL**
Thigh circumference*(cm)	30 - 36	36 - 44	44 - 50	50 - 57	57 - 64	64 - 72	72 - 78

* measured 15 cm above patella centre

** available only for wrap-around version

Length: 31 cm

MEYRA
medical CE

Stabilizing knee brace

Type: without regulation of flexion movement

Catalogue number: QKN01021, QKN01031,
QKN01041, QKN01051

PRODUCT CHARACTERISTICS:

- made of polyamide knitted fabric laminated with a layer of perforated neoprene,
- polycentric stays provide lateral stabilization,
- special construction allows to adjust the stabilization strength to individual needs,
- patellar opening and opening for popliteal fossa which enable unhampered knee flexion.

USE:

Intense, chronic pain syndromes of relatively low intensity and of various etiology, postoperative stabilization during final phase of functional treatment, injuries and dysfunctions which do not require joint immobilization, support for the ligaments (apparatus) of knee joint.

Size	S	M	L	XL
Thigh circumference* (cm)	41 - 44	44 - 47	47 - 50	50 - 53

* In order to choose a proper size, measure the thigh circumference 15 cm above the center of patella

Length: 31 cm



Qmed® Silver Line CE

Stabilizing knee brace

Type: with regulation of flexion movement

Catalogue number: QKN04020, QKN04030,
QKN04040, QKN04050

PRODUCT CHARACTERISTICS:

- made of polyamide knitted fabric laminated with neoprene,
- stays with polycentric hinges provide lateral stabilization and allow to control flexion/extension movement,
- patellar opening and straps enable optimal adjustment and regulation of stabilization strength,
- patellar opening and opening for popliteal fossa which enable unhampered knee flexion,
- flexion/extension stops 15°, 30°, 60°, 90°.

USE:

Intense, chronic pain syndromes of relatively low intensity and of various etiology, functional treatment of knee joint injuries which require control of angular range of motion, injuries and dysfunctions which do not require joint immobilization, support for the ligament apparatus of knee joint.

Size	S	M	L	XL
Thigh circumference* (cm)	41 - 44	44 - 47	47 - 50	50 - 53

* In order to choose a proper size, measure the thigh circumference 15 cm above the center of patella

Length: 31 cm



Qmed® Silver Line CE

S-move open

Open knee brace

Catalogue number: DRQI01



PRODUCT CHARACTERISTICS:

- made of advanced, multi-directionally flexible, breathable fabric AirX® with Clariant's moisture managing finish,
- quick drying, maintains its original features after washing,
- main body made of knitted fabric adjusts to anatomical shape of the limb and patella,
- open construction enables easy and quick putting on as well as width adjustment,
- two aluminium stays enable lateral stabilization,
- four axis, patented hinge,
- regulated extension (0°, 10°, 20°, 30°) and flexion (45°, 60°, 90°),
- additional stabilizing straps of opposite courses enable proper position of stays and prevent the brace from uncontrolled movement and rotation.

57

USE:

Overstrain and micro-injuries, subacute and chronic pain syndromes of different aetiologies, degenerative changes, knee laxity, further treatment of knee lateral ligaments.

Size	1	2	3	4	5
Shin circumference* (cm)	36 - 39	40 - 43	44 - 48	49 - 54	55 - 60

* In order to choose a proper size, measure the circumference 15 cm above the center of patella

Length: 34,5 cm

Qmed® CE

S-move

Stabilizing knee joint orthosis

Catalogue number: DRQI0D



PRODUCT CHARACTERISTICS:

- made of advanced, multi-directionally flexible, breathable fabric AirX®,
- polycentric stays for lateral stabilization,
- hinges for double limits to the range of extension (0°, 10°, 20°, 30°) and flexion (45°, 60°, 90°),
- two additional stabilising straps allow to adjust compression,
- shaped patellar opening provides lateral stabilisation and relief of the femoropatellar joint.

USE:

Stabilisation after surgery in the final phase of functional treatment, additional treatment for knee twists, stabilisation in the case of degenerative changes, acute, subacute and chronic pain syndromes of different aetiologies. Injury and dysfunction without joint immobilisation, knee ligament system supports in overstraining and instability, relief for ligamentous-muscular system in overloads, patella stabilisation.

Size	1	2	3	4	5
Knee circumference (cm)*	36 - 39	40 - 43	44 - 48	49 - 54	55 - 60

* In order to choose a proper size, measure the circumference 15 cm above the center of patella

Length: 35,5 cm

Qmed® CE





Motive Plus

Open knee brace

Catalogue number: DRQ105

60



PRODUCT CHARACTERISTICS:

- open knee brace made of elastic, breathable spacer fabric AirX®,
- equipped with patellar opening,
- covered popliteal fossa,
- stays with polycentric hinge provide lateral stabilization,
- transverse straps provide additional stabilization,
- colour: black.

USE:

Repetitive, subacute pain related to strain injury and micro-injury of muscles and ligaments, gonarthrosis and knee joint instability, stabilization after injury during the final phase of functional treatment, injury and dysfunction which do not require knee immobilization (lateral ligament injuries).

Size	1	2	3	4	5
Thigh circumference* (cm)	36 - 39	40 - 43	44 - 48	49 - 54	55 - 60

* In order to choose a proper size, measure the thigh circumference 15 cm above the center of patella

Length: 40 cm

Qmed® CE

Reflex Long

Knee ROM brace

Catalogue number: Wrap-around version: MQB3N

Closed version: MQB3O



PRODUCT CHARACTERISTICS:

- made of breathable spacer fabric which carries humidity away and provides optimal comfort of use,
- lateral splints may be adjusted to anatomical shape of the limb which allows to control range of motion and improves comfort of use,
- removable straps may be applied on right and left side and their marking facilitates proper fastening,
- delicate, crease-resistant fabric of specially chosen tension protects against irritation and compression of popliteal fossa,
- closed version equipped with patella stabilizer ring which location and length can be adjusted according to patient's need,
- colour: black.

USE:

Mild knee instability of various etiology.

Size	XS	S	M	L	XL	XXL	XXXL**
Thigh circumference*(cm)	30 - 36	36 - 44	44 - 50	50 - 57	57 - 64	64 - 72	72 - 78

* measured 15 cm above patella centre

** available only for wrap-around version

Length: 40 cm

MEYRA *medical* CE

GENUmeyra IMS

Knee immobilisation splint

Catalogue number: MY2712



PRODUCT CHARACTERISTICS:

- padded knee splint with soft edge,
- anatomically formed reinforcement,
- easy and steady fastening by means of hook and loop straps,
- removable pad for the hollow of the knee,
- left and right version,
- colour: navy blue.

USE:

Postoperative immobilization of the knee joint, ligament and kneecap injuries.

Size	0	1	2
Length (cm)	40	48	58

MEYRA
medical CE

Shin and thigh brace

Catalogue number: QKM010



PRODUCT CHARACTERISTICS:

- made of velour polyamide knitted fabric,
- soft padding provides optimal comfort of use,
- three stays (one rear splint and two symmetrical, lateral stays) ensure knee joint immobilization,
- thanks to additional straps and possibility of changing the position of lateral stays, the brace optimally adjusts to body shape,
- patellar opening.

USE:

Conditions which require knee joint immobilization after surgeries or injuries.

Size	S	M	L	XL	XXL
Thigh circumference* (cm)	41 - 44	44 - 47	47 - 50	50 - 53	53 <

* In order to choose a proper size, measure the thigh circumference 15 cm above the center of patella

Length: 56 cm

Qmed[®] Silver Line CE

Post-operative knee brace with hinge

Catalogue number: QKM06090

PRODUCT CHARACTERISTICS:

- two splints provide lateral stabilization,
- light, open construction protects against skin abrasions,
- four self-adhesive stabilizations straps,
- easy to use hinge to control flexion (0°-120°) and extension (0°-90°) movements every 10°,
- one, universal model for left and right limb.

USE:

Conditions which require limitation of flexion/extension movement after surgery or injury.

One, universal size
Length from 49,5 cm to 62 cm



Qmed® Silver Line CE

ROM Knee brace covering shin and thigh

Catalogue number: DRQ135

PRODUCT CHARACTERISTICS:

- splints for lateral stabilization,
- light, open construction prevents skin abrasion and rash,
- four stabilizing straps with hook and loop closure,
- easy to operate hinge, adjustable range of motion (flexion: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°, 105°, 120°, extension: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°),
- quick lock 0°, 15°, 30°, 45°,
- universal model for right and left limb.

USE:

Conditions which require immobilization or control of knee joint range of motion after surgeries or injuries.

One, universal size
Adjustable length from 56 cm to 72 cm



Qmed® CE

Rush

Post-operative knee brace

Catalogue number: DRQ10N

PRODUCT CHARACTERISTICS:

- two splints provide lateral stabilization,
- five self-adhesive stabilizing straps (including infrapatellar strap), each independently adjusted,
- easy to operate monocentric hinge which controls range of motion every 10°, flexion: 0°- 120°, extension: -10°- 90°,
- quick lock: 0°, 15°, 30°, 45°,
- flat hinge,
- intuitive length adjustment,
- front panel enables quick and easy adjustment of brace circumference,
- extra option: cold pack which can be attached to the hinge,
- soft padding protects the leg and ensures its optimal position,
- one, universal model for left and right leg,
- weight of the brace does not exceed 800 grams.

USE:

Conditions which require immobilization or control of knee range of movement after surgeries or injuries. Treatment and further treatment of knee joint and soft tissue injuries.

One, universal size

Adjustable length: 48 – 61 cm (thigh: 22 – 28 cm, shin: 26 – 33 cm)



Qmed® CE

Compact Pro

Post-operative knee brace

Catalogue number: OV51600

64



PRODUCT CHARACTERISTICS:

- two splints provide lateral stabilization,
- four self-adhesive stabilizing straps, each independently adjusted,
- easy to operate monocentric hinge which controls range of motion: flexion: 0° - 120°, extension: 0° - 90°,
- quick lock: 0°, 15°, 30°, 45°,
- flat hinge,
- intuitive length adjustment,
- front panel enables quick and easy adjustment of brace circumference,
- soft padding protects the leg and ensures optimal position,
- one, universal model for left and right leg.

USE:

Conditions which require immobilization or control of knee range of movement after surgeries or injuries. Treatment and further treatment of knee joint and soft tissue injuries.

One, universal size

Adjustable length from 49 cm to 66 cm

MEYRA
medical CE

Follow

Knee brace with regulated range of motion

Catalogue number: DRQ10E



PRODUCT CHARACTERISTICS:

- anatomically profiled, light construction,
- frame made of high resistance aluminium commonly used in aircraft industry,
- patented, four-axis hinge which reflects physiological knee movement (rolling and sliding),
- four-point support system with 0° immobilization,
- regulation of extension 10°, 20°, 30° and flexion 10°, 20°, 30°, 45°, 60°, 90°, 110°,
- movable attachments for the straps,
- made of stainless materials, resistant to water, salt, low temperature.

USE:

Decreased post-traumatic capacity of the muscle and ligament structure that stabilize the knee joint, condition after reconstructing the ACL, PCL, MCL and LCL ligaments, ACL and PCL ligament failure, weakened strength of the quadriceps muscle, regaining sports activity after knee joint injuries.

Size	XS	S	M	L	XL
Thigh circumference* (cm)	33-39	39-47	47-53	53-59	59-67
Knee circumference (cm)	30-33	33-35	35-38	38-43	43-48
Calf circumference** (cm)	25-30	30-35	35-40	40-47	47-50

* measured 15 cm above patella

** measured 15 cm below patella

Qmed[®] CE

Game Changer

OA brace

Catalogue number: left - OV 49001
right - OV 49002



PRODUCT CHARACTERISTICS:

- light, perforated, elastic frame,
- pads made of modern fabric allow to position the brace properly (anti-migration system with silicone pads),
- one time adjustable system of stabilizing straps, easy to use buckles,
- repetitive, easy to measure system of corrective fore adjustment,
- flat and light hinge with hyperextension blockade,
- two versions available: left and right.

USE:

Uni-compartment knee osteoarthritis, further postoperative treatment, rehabilitation process which require offloading of the degenerative compartment, correction and stabilization of knee joint, useful for patients of limited and medium activity level.

One, universal size
Length: 35 cm



OvationMedical®



First

Ligament knee brace

Catalogue number: DRQ15A



PRODUCT CHARACTERISTICS:

- anatomically profiled, light aluminium frame,
- polycentric hinge provides flexion control (0°, 15°, 30°, 45°, 60°, 75°, 90°), extension control (0°, 10°, 20°, 30°, 40°) and 0° immobilization,
- stabilization straps keep the brace in place and ensure effective treatment,
- breathable, condylar pads placed below the hinge protect against skin abrasions and provide maximum comfort,
- two extra padding elements for change,
- two versions available: for right and left leg.

USE:

Stabilization after ACL, PCL, MCL, LCL, knee joint protection against recurrent injuries while returning to sport activities, follow-up treatment after reconstruction of ligaments in knee joint.

Size	XS	S	M	L
Thigh circumference* (cm)	33-39	39 - 47	47 - 53	53 - 60

* In order to choose a proper size, measure the thigh circumference 15 cm above the center of patella

Length: 35 cm size S, XS
36 cm size M, L

Qmed®





Grand
Follow

OVAL
ACTIVE



KAFO type brace

Catalogue number: DRQI5D

68



PRODUCT CHARACTERISTICS:

- four self-adhesive stabilizing straps,
- light, open construction prevents abrasions and sore skin,
- ankle joint stabilizer (AFO) made from aluminium and EVA,
- ankle joint stabilizer (AFO) of adjustable height,
- two aluminium splints for lateral stabilization,
- easy to use, monocentric hinge: flexion: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°, 105°, 120°, extension: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°,
- quick lock: 0°, 15°, 30°, 45°.

USE:

Conditions which require complete immobilization of lower limb intra-articular fractures, fractures of tibia, fibula and femur.


SIZE CHART:

Length without ankle joint stabilizer (AFO): 56 - 72 cm (A)

Length of ankle joint stabilizer (AFO): 13,5 - 16 cm (B)

Length with ankle joint stabilizer (AFO): 69,5 - 88 cm (C)

Size	S	M	L
Foot pad length (cm) (D)	16	17,5	19
Foot length (cm)	<22	22-25	>25



Qmed[®] CE

Hip joint orthosis

Catalogue number: DRQI0B



PRODUCT CHARACTERISTICS:

- consists of elastic straps which cover hips and thigh band made of polyamide knitted fabric, finished with hook and loop fastener,
- hip pad made of the polyurethane foam laminated with velour knitted fabric,
- pentapro foil edged with the velour knitted fabric works as additional, removable stabilization element,
- stabilizes hip joint, enables controlled motion and provides compression of soft tissues,
- left and right version,
- colour: black.

USE:

Help in treatment and rehabilitation process in case of injuries, illnesses and hip joint dysfunctions, post-operative hip compression after anastomoses of intertrochanteric fractures, intertrochanteric osteotomies and the endoprosthetics, overstrain changes of the hip joint, degenerative changes of the hip joint, pain syndromes of various etiology.

Size	S	M	L	XL
Hip circumference (cm)	75 - 90	90 - 105	105 - 120	120 - 140

Qmed[®] CE

Ankle support

Catalogue number: PHAI011



PRODUCT CHARACTERISTICS:

- multifunctional ankle support,
- protects and supports knee joint,
- provides optimal joint compression,
- ideal for practicing sport,
- helpful in rehabilitation process.

USE:

Strain injury, reinforcement of ankle joint, protection against injuries.

Size	S	M	L	XL
Ankle circumference (cm)	19 - 21	22 - 24	25 - 27	28 - 30

Qmed+pharma CE

ACHILLOmedical

Achilles tendon bandage with gel pads

Catalogue number: MY27212




PRODUCT CHARACTERISTICS:

- knitted fabric for optimal comfort and better adjustment to body shape,
- lateral silicone pads relieve Achilles tendon,
- silicone heel wedges for tension reduction,
- left and right version,
- colour: silver-blue.

USE:

Ruptures and injuries of the Achilles tendon, support of the Achilles tendon after ruptures, postoperative.

Size	1	2	3	4	5
Ankle circumference (cm)	17 - 19	19 - 21	21 - 23	23 - 25	25 - 27



MEYRA medical CE





MALLEOmedical

Ankle joint bandage with gel pads

Catalogue number: MY2722

PRODUCT CHARACTERISTICS:

- knitted fabric for optimal comfort and better adjustment to body shape,
- lateral silicone pads for better compression pressure distribution,
- elastic edge without constrictions,
- left and right version,
- colour: silver-blue.


USE:

Ligament injuries or weakness, arthrosis and distorsions, tendomyopathies, irritations.

Size	1	2	3	4	5	6
Ankle circumference (cm)	17 – 19	19 – 21	21 – 23	23 – 25	25 – 27	27 – 29



MEYRA
medical



AFO - Soft

Stabilizing ankle joint orthosis

Catalogue number: DRQ10H

PRODUCT CHARACTERISTICS:

- open type orthosis, made of Air-X spacer fabric of very good air and moisture permeability,
- "comfort areas" placed on vulnerable spots in order to protect against abrasions and provide better stabilization,
- hook-and-loop closure,
- additional, flexible, stabilizing strap in order to control compression and provide proper foot position,
- universal version for left and right leg.

USE:

Effectively supports treatment and rehabilitation of injuries, ankle joint dysfunctions, provides additional help to the lateral ankle ligaments, as well as medial orientation of the foot.

Size	S	M	L
Leg circumference above the ankle (cm)	19 - 22	23 - 25	26 - 28



Qmed® CE

Stabilizing ankle brace (pneumatic-gel)

Catalogue number: QAR03091



PRODUCT CHARACTERISTICS:

- plantar part of regulated length allows to adjust the circumference and ensures unhampered dorsiflexion/plantarflexion movement,
- internal part of the brace consists of two elements: upper – filled with air and bottom part filled with gel (cooling/warming).

USE:

Post-traumatic conditions of ankle joint (sprain, dislocation) during functional treatment, instability of ankle joint, stabilization after surgeries, ankle joint protection from recurrent injuries while re-undertaking occupational or sport activities after experiencing injury.

One, universal size

73

Qmed® Silver Line CE

TALUmedical

Stabilising ankle joint orthosis

Catalogue number: MY27213



PRODUCT CHARACTERISTICS:


- preformed stable plastic shells for a better fit,
- two hook and loop fasteners for a steady fixation,
- adjustable foot rest,
- foam pad for an optimal adjustment.

USE:

Torn ligaments, sprains and distorsions of the upper ankle joint, postoperative after treatment.

One, universal size

MEYRA
medical CE



AeroStab

Stabilizing ankle brace

Catalogue number: DRQI0P



PRODUCT CHARACTERISTICS:

- made of spacer fabric of very good air and moisture permeability,
- intuitive system of lace fastening adjusts compression to patient's individual needs,
- anatomically shaped stays, additionally reinforced with thermoplastic material, provide optimal stabilization,
- soft padding protects soft tissues,
- one, universal model for left and right leg,
- colour: black.

USE:

Sprain of tibio-tarsal joint (I° i II°), protection against excessive inversion and eversion, habitual sprain of tibio-tarsal joint, support in recovery process, prevention of injuries.

Size	S	M	L
Ankle circumference (measured over the heel)	27-30	30-33	33-36

Qmed® CE

AFO - Walker 3/4

Step-ankle orthosis

Catalogue number: DRQI3H



PRODUCT CHARACTERISTICS:

- short step-ankle orthosis of light, rigid and durable structure,
- covered hinge enables limitation of plantar and dorsal flexion (0°-70°) every 10°,
- removable, soft lining,
- fastened by self-adhesive band,
- one, universal version for left and right leg.

USE:

Protection and stabilization of foot and ankle after injuries, necessity to immobilize ankle joint in the chosen therapeutic range, postoperative, preventive treatment.

Size	S	M	L	XL
Footwear size	34 - 39	39 - 44	44 - 47	47 - 50
Orthosis height (cm)	25,4	25,4	25,4	25,4

Qmed® CE

Foot and shin brace

Catalogue number: QAM04010, QAM04020,
QAM04030, QAM04040, QAM04050

76



PRODUCT CHARACTERISTICS:

- equipped with velour knitted-fabric padded with foam,
- five additional hook-and-loop closures which allow to adjust compression to individual needs,
- two stays provide lateral stabilization,
- hinge allows to regulate the dorsiflexion / plantarflexion (0° - 45°) every 15°,
- two versions: left and right,
- anatomically moulded outsole enables stabile, similar to physiological walk.

USE:

Fractures/cracks of metatarsus bones also during post-operative treatment, cracks of heel/tarsus/ great toe bones, surgeries of great toe joint, soft tissue injury, injuries of calcaneal tendon, sprain/dislocation of ankle joint, protection and post-traumatic stabilization of foot and ankle joint area, immobilization of ankle joint according to chosen therapeutic range, post-operative treatment.

Size	XS	S	M	L	XL
Foot length (cm)	< 21	21 - 24	24 - 27	27 - 30	30 - 33

Length: size XS - 24 cm, size S/M - 43 cm, size L/XL - 45 cm

Qmed® Silver Line CE

Foot and shin brace (pneumatic)

Catalogue number: QAM04011, QAM04021, QAM04031,
QAM04041, QAM04051



PRODUCT CHARACTERISTICS:

- soft, pneumatic padding allows to control the soft tissue compression which prevents oedma and provides optimal adjustment to body shape,
- fastened with hook-and-loop closure,
- hinge allows to regulate the dorsiflexion / plantarflexion (0° - 45°) every 15°,
- anatomically moulded outsole enables stabile, similar to regular walk.

USE:

Stabilization and counteraction against oedma in case of soft tissue injuries of ankle joint/distal part of shin, a second/third degree ankle sprain, follow-up treatment in case of stabile fractures of ankle/instep, immobilization in required therapeutical position after surgeries, surgeries of great toe joint, soft tissue injury, injuries of calcaneal tendon, protection and post-traumatic stabilization of foot and ankle joint area, immobilization of ankle joint according to chosen therapeutic range.

Size	XS	S	M	L	XL
Foot length (cm)	< 21	21 - 24	24 - 27	27 - 30	30 - 33

Length: size XS - 24 cm, size S/M - 43 cm, size L/XL - 45 cm

Qmed® Silver Line CE

Pneumatic Walking Brace

Foot and shin brace

Catalogue number: OV1000



PRODUCT CHARACTERISTICS:

- frame made from fiberglass reinforced with nylon,
- liner made from open-cell foam with brushed nylon,
- latex-free,
- edgeless sole ensures 3-Dimensional gait and easy pivoting,
- colour: black and blue,
- one, universal model for left and right leg.

USE:

Soft tissue injuries of ankle, foot and shin, ankle sprains (I, II degree), further treatment of stable fractures of metatarsus/ankle.

CONTRAINDICATIONS:

Disorder of sensation, skin irritations, open wounds, numbness, pain.

Length: 43 cm

Size	Shoe size
XS	34-36
S	37-39
M	40-44
L	45-47



OvationMedical®



Pneumatic Short Walking Brace

Foot and shin brace (short)

Catalogue number: OV1010



PRODUCT CHARACTERISTICS:

- frame made from fiberglass reinforced with nylon,
- liner made from open-cell foam with brushed nylon,
- latex-free,
- edgeless sole ensures 3-Dimensional gait and easy pivoting,
- colour: black and blue,
- one, universal model for left and right leg.

USE:

Soft tissue injuries of ankle, foot and shin, further treatment of stable fractures of metatarsus/ankle.

CONTRAINDICATIONS:

Disorder of sensation, skin irritations, open wounds, numbness, pain.

Length: 32 cm

Size	Shoe size
XS	34-36
S	37-39
M	40-44
L	45-47



OvationMedical®





Air Walking Boot

Ankle foot orthosis

Catalogue number: DRQ15F



PRODUCT CHARACTERISTICS:

- rigid, perforated cover provides effective, external stabilization,
- rigid cover ensures air circulation,
- profiled outsole enables stabile and physiological walk,
- inner lining made of soft, breathable fabrics,
- two air chambers allow to adjust the soft tissue compression and stabilization force,
- adjustable cover for toes protection,
- universal right and left limb supply.

USE:

Air Walking Boot provides stabilization and protection of foot and ankle joint area. Follow-up treatment in case of stabile fractures of ankle/instep.

CONTRAINDICATIONS:

Open, unstable fractures; difficult to heal wounds, sores and other skin changes.

SIZE:

Short brace	S	M	L	XL
Height (cm)	19,5	21,5	22,5	25
Insole size	36-38	38-40	40-43	43-46
Weight (kg)	0,8	0,93	1,05	1,40

Qmed® CE

Air Walking Boot

Ankle foot orthosis

Catalogue number: DRQ15E



PRODUCT CHARACTERISTICS:

- rigid, perforated cover provides effective, external stabilization,
- rigid cover ensures air circulation,
- profiled outsole enables stabile and physiological walk,
- inner lining made of soft, breathable fabrics,
- two air chambers allow to adjust the soft tissue compression and stabilization force,
- adjustable cover for toes protection,
- universal right and left limb supply.

USE:

Air Walking Boot provides stabilization and protection of foot and ankle joint area. Follow-up treatment in case of stabile fractures of ankle/instep.

CONTRAINDICATIONS:

Open, unstable fractures; difficult to heal wounds, sores and other skin changes.

SIZE:

High brace	S	M	L	XL
Height (cm)	33	37	39	41
Insole size	36-38	38-40	40-43	43-46
Weight (kg)	1,00	1,16	1,32	1,76

Qmed® CE

Power Walking Boot

Tibial Ankle Foot Orthosis

Catalogue number: Q62011

80



PRODUCT CHARACTERISTICS:

- rigid, perforated construction which provides effective stabilization and enables free air circulation,
- profiled sole enables stable, similar to regular walk,
- inner lining made of soft, breathable fabric,
- additional toe protection,
- internal air spaces allow to adjust soft tissue compression and choose appropriate stabilization,
- universal for right and left limb.

USE:

Protection against oedema, stabilization in case of: soft tissue injuries (ankle area, foot distal part of tibia), sprain of the ankle, metatarsus, tarsus (second and third degree). Further treatment of stable metatarsus and ankle fractures, post-operative treatment, inflammation of Achilles tendon, extensor tendon, foot flexor.

Size	S	M	L	XL
Footwear size	36 - 39	38 - 42	42 - 46	> 46

Length: S - 42 cm, M - 46 cm, L/XL - 48 cm

Qmed[®] Silver Line CE

Offloading insoles

Catalogue number: DRQR5A1



PRODUCT CHARACTERISTICS:

- made of EVA foam covered with fleece,
- five sizes in one set,
- two versions available: left and right,
- colour: white.

USE:

Treatment and rehabilitation of Achilles tendon injuries.

Qmed[®] CE





Scamper

Pneumatic boot with compression control

Catalogue number: DRQI3T



PRODUCT CHARACTERISTICS:

- rigid, perforated shell provides stabilization and allows free air circulation,
- soft padding with air cells for soft tissue compression,
- unique, patented control system to adjust and regulate compression,
- prevents uncontrolled, too strong compression and protects against tissue trophism,
- anatomically moulded outsole enables stable, physiological walk,
- universal model for left and right leg.

USE:

Stabilization and counteraction against oedma in case of soft tissue injuries of ankle; second or third degree ankle sprain; follow-up treatment in case of stabile fractures of instep or ankle; immobilization after surgeries; surgeries of great toe joint; soft tissue injury; injuries of calcaneal tendon; protection and post-traumatic stabilization of foot and ankle joint.

Size	S	M	L	XL
Shoe size (cm)	36-37	38-42	42-46	> 46

Qmed® CE



AFO - Fit

Foot-drop orthosis

Catalogue number: DRQI1A



PRODUCT CHARACTERISTICS:

- light,
- made of polypropylene,
- additional band to stabilize ankle joint,
- free heel and calf,
- highly hygienic,
- two versions: for left and right leg.

USE:

Prevention of foot-drop in the case of peroneal nerve palsy and neurosurgical diseases.

Size	S	M	L	XL
Footwear size	35 - 37	37 - 39	39 - 41	41 - 44

Qmed® CE

Hybrid Night Splint

Foot Drop Brace

Catalogue number: OV3001

84



PRODUCT CHARACTERISTICS:

- frame made from polypropylene,
- liner made from an opencell foam with brushed nylon,
- latex free,
- fully adjustable tension straps and articulating hinge design provide optimal support and comfort of use,
- colour: black and blue,
- one, universal version for left and right leg.

USE:

Inflammation of plantar aponeurosis, peroneal nerve injury, post-operative conditions of shin and foot.

CONTRAINDICATIONS:

Varicose veins, skin irritation, swollen shin and foot, diabetic foot, congenital disorders, spasticity, congenital foot deformities.

Size	Shoe size (Men)	Shoe size (women)
S-M	39-42,5	37,5-42,5
L-XL	43-49	43-48



Ovation Medical®



AFO - Night

Foot - drop orthosis

Catalogue number: DRQ13A



PRODUCT CHARACTERISTICS:

- light structure made of polypropylene,
- system of stabilizing stripes with soft lining enables correction of foot position to the shin,
- open type orthosis with spacious interior,
- possible to use with a cast on foot and shin.


USE:

Correction of feet position to the shin from medial position to dorsal bend of the foot, prevention of contractures in bedridden patients with hemiplegia, multiple sclerosis, cerebral palsy.

Size	S	M	L
Length (cm)	23	24,5	26

Qmed®





WWW.MDH.PL

WWW.QMEDINFO.PL

MDH SP. Z O.O.
UL. KS. W. TYMIENIECKIEGO 22/24
90-349 ŁÓDŹ
POLAND

MARKETING DEPARTMENT:
TEL. 42/636 54 82
E-MAIL: MARKETING@MDH.PL

SALES DEPARTMENT:
TEL. 42 636 81 44
FAX 42 299 69 24
E-MAIL: EXPORT@MDH.PL